

ophthalmic

Contents

Reusable

Scissors	1 - 10
Surgical Knives and Blades	11
Diamond Knives	12 - 14
Hooks, Probes and Chops	15 - 27
Irrigation / Aspiration	28 - 30
Specula	31 - 35
Retractors and Rings	36 - 37
Needleholders	38 - 45
DCR Instruments	46 - 51
CDCR Tubes - Lester Jones, Cox System	52 - 53
Depressors, Punches and Calipers	54 - 55
Pressure Reducers and Burrs	56
Bipolar Forceps and Pencils	57 - 58
Bipolar Cables	59
Forceps	60 - 80
Lid Surgery	81 - 84
Corneal Markers and Global Fixation	85 - 86
Corneal Trephines and Punches	87
Enucleation and Orbit	88
Vitreoretinal Instruments	89
Volk Lenses	90 - 92
Microwash Trays and Pillows	93
Sterilization Trays	94
Surgeons Stools and Q-Rinse	95

Single Use

Corneal Trephines and Punches	98 - 99
DCR	100
Oculoplastic Oncology	101
Bipolar	102
Vitreoretinal	103 - 109
Volk Lenses	110
Punctal Occlusion	111 - 112
Blades and Knives	113 - 114
Mani Knives	115 - 118
Cellulose and PVA Products	119 - 120
Miscellaneous	121 - 122
Sterimedix Cannula	123 - 142
Sterimedix Coaxial I/A	133 - 134
Sterimedix Bimanual I/A	135 - 136
Drapes and Fluid Collection	143

Repair Service

We offer a cost effective repair service for any type and brand of ophthalmic surgical instrument, and for most types of phacoemulsification handpieces.

Instructions for use.

Reprocessing instructions.

Product instructions for use, decontamination instructions, plus various product statements and company certificates are available www.altomed.com.

Service

Altomed Area Managers, based throughout the UK are available to visit you at your location to demonstrate products, offer technical advice or help with administrative issues.

Microwash Custom Trays

We are able to offer a totally bespoke instrument tray solution to help you keep your delicate instruments safe and secure during decontamination processes.

For information about Microwash see page 93, or go to www.microwash.net.

Visit www.altomed.com



Email: sales@altomed.com Web: www.altomed.com

Tel: +44 (0)191 519 0111 Fax: +44 (0)191 519 0283

Altomed Ltd, 2 Witney Way, Boldon Business Park, Tyne & Wear, NE35 9PE. England

Reg. No: 01243526

Scissors

Eye Scissors

A1060

Pointed tips,
straight
90mm



A1065

Pointed tips,
curved
90mm



A1070

Rounded tips,
straight
90mm



A1075

Rounded tips,
curved
90mm



Iris Scissors

A1100

Straight
90mm



A1105

Curved
90mm



Scissors

Iris Scissors

A1110

Straight
100mm



A1115

Curved
100mm



A1120

Straight
115mm



A1125

Curved
115mm



Strabismus Scissors

A1130

Straight
115mm



A1135

Curved
115mm



Scissors

Strabismus Scissors

A1140

Straight
90mm



A1145

Curved
90mm



Stevens Tenotomy Scissors

A1150

Pointed tips,
straight
115mm



A1155

Pointed tips,
curved
115mm



A1160

Rounded tips,
straight
115mm



Scissors

Stevens Tenotomy Scissors

A1165

Rounded tips,
curved
115mm



Pooley Conjunctival Scissors

A1170

Pooley conjunctival
scissors
90mm



Werbs Scissors

A1180

Werbs scissors
115mm



Enucleation Scissors

A1185

Enucleation scissors
125mm



Scissors

Bonn Scissors

A1190

Straight
90mm



A1195

Curved
90mm

Stitch Scissors

A1235A

Stitch scissors
90mm



Dressing Scissors

A1240

Rounded tips
150mm



A1244

Pointed tips
150mm



Scissors

Conjunctival Scissors

A1250

Rounded tips,
curved
115mm



A1251

Rounded tips,
straight
115mm



A1252


Pointed tips,
straight
115mm



Williamson Noble Scissors


A1300

Pointed tips,
straight
120mm



A1310

Rounded tips,
straight
120mm



A1312

Pointed tips,
curved
120mm



A1314

Rounded tips,
curved
120mm



Scissors

De Wecker Scissors


A1315

Pointed tips
110mm



A1325

Rounded tips
110mm



Castroviejo Scissors

A1350

Rounded tips
100mm



A1355

Pointed tips
100mm



Westcott Tenotomy Scissors


A1360

Pointed tips
115mm



A1365

Rounded tips
115mm



Scissors

Vannas Scissors

A1380

Short blades,
straight
80mm



A1385

Short blades,
curved
80mm




A1395

Short blades,
angled
80mm



A1400

Long blades,
straight
85mm




A1405

Long blades,
curved
85mm



A1415

Long blades,
angled
85mm



Scissors

Capsule Scissors


A1435

Ong type,
angled
80mm



A1440

Round handle,
curved
110mm



A1445


Flat handle,
curved
110mm



Castroviejo Corneal Scissors


A1472

6mm blades,
left
115mm





A1473

6mm blades,
right
115mm



A1490

6mm blades,
left
100mm



A1495

6mm blades,
right
100mm



Scissors

Katzin Corneal Scissors

A1510

9mm blades, left
100mm



A1515

9mm blades, right
100mm



Castroviejo Micro Scissors

A1520

Pointed tips,
straight
100mm



A1525

Rounded tips,
curved
100mm



A1530

Rounded tips,
strong curve
100mm



Surgical Knives and Blades

A1650

Daya lamellar
dissector



A1648

Barkan goniotomy
knife



A1660

Tooke knife,
straight



A1662

Tooke knife,
angled



A1696

Troutman razor
fragment holder,
straight
120mm



A1702

Scalpel handle, No.3



A1725


Metal handle, for
Single Use blades



Free hand diamond knives with straight handle for cataract surgery


Only a small selection of Meyco Diamond knives are shown here

Please contact customer services for details of the full range



ME-100

45° Single edged diamond knife




ME-100

B	L	D	A	S
1.50	3.50	0.20	45°	1.50

ME-106

30° Single edged diamond knife




ME-106

B	L	D	A	S
1.50	3.50	0.20	30°	1.70

ME-120

Double lancet diamond knife




ME-120

B	L	D	A	S
1.50	3.50	0.20	60°	1.30

ME-130

Tri-facet diamond knife



ME-130

B	L	D	S
1.50	3.50	0.20	1.30°



Diamond Knives

ME-231

2.5mm Phaco knife,
angled

ME-232

2.65mm Phaco knife,
angled

ME-233

2.8mm Phaco knife,
angled

ME-235

3.0mm Phaco knife,
angled

ME-237


3.2mm Phaco knife,
angled

Phaco knives with angled handle for
"clear cornea" technique



Double lancet 4 mm
blade with sharp sides

	B	L	D	A	S
ME-231	2.50	4.00	0.20	90°	1.25
ME-232	2.65	4.00	0.20	90°	1.30
ME-233	2.80	4.00	0.20	90°	1.40
ME-235	3.00	4.00	0.20	90°	1.50
ME-237	3.20	4.00	0.20	90°	1.60



Knife for scleral tunnel incision


ME-285

2mm Scleral tunnel knife



ME-285 Tunnel Blade


B	L	D
2.00	4.00	0.20



Limbal relaxing incision knife


ME-701

1mm L.R.I. knife, 0.65mm
cutting depth



ME-701


Cutting depth	B	B1	L	D	A
0.65mm	1.00	0.40	2.00	0.20	40°



Diamond Knives



Diamond knives with micrometer for refractive surgery




ME-420 Type RK-T

B	L	D
1.00	3.50	0.20

ME-420

Thornton RK-T knife
(Inc. measuring gauge)




ME-425 Type RK-TC

B	L	D	A	S
0.80	3.50	0.10	15°	0.20

ME-425

Thornton RK-TC knife
(Inc. measuring gauge)



ME-451 Type RK-M


B	L	D	A	S
1.00	3.50	0.10	45°	1.00

ME-451

45° RK-M knife
(Inc. measuring gauge)

Calibration gauge for checking of the depth of incision setting.

We recommend that the cutting depth is checked before every procedure. The cutting depth can be measured using an ME-409 (for DK-Blades double edged) or ME-459 (for M-Blades single edged).



Hooks, Probes and Chops

A2012

Anel lacrimal probes,
set of three sizes



A2020

Bowman lacrimal
probe, 0.9mm x
1.0mm, 00/0



A2022

Bowman lacrimal
probe, 1.1mm x
1.2mm, 1/2



A2024

Bowman lacrimal
probe, 1.3mm x
1.4mm, 3/4



A2026

Bowman lacrimal
probe, 1.5mm x
1.6mm, 5/6



A2018

Bowman lacrimal
probes, set of four
sizes
(As shown above)

Hooks, Probes and Chops

A2040

Liebreich lacrimal probe, 0.9 x 1.0mm, 00/0



A2042

Liebreich lacrimal probe, 1.1 x 1.2mm, 1/2



A2044

Liebreich lacrimal probe, 1.3 x 1.4mm, 3/4



A2046

Liebreich lacrimal probe, 1.5 x 1.6mm, 5/6



A2048

Liebreich lacrimal probe, extra small, 000/0000



A2038

Liebreich lacrimal probes, set of five (as shown above)

Hooks, Probes and Chops

A2050

Pigtail probe
140mm



A2080

Nettleship dilator



A2087

Hays punctum finder



A2091

Infant
dilator



A2094

Spud and Gouge
needle (with caps)



Hooks, Probes and Chops


A2110

Fosters D/E lacrimal dilator and finder (with cap)



A2122

Fosters lacrimal probes, set of three, 0.6, 0.8 and 1.0mm



A2135

Kimura spatula 150mm



A2142

Castroviejo spatula, small 10mm



A2144

Castroviejo spatula, large 15mm



A2145

Castroviejo double ended spatula, 10mm / 15mm



Images not to scale

Hooks, Probes and Chops

A2152

Sinclair
spatula



A2154

Nettleship iris
repositor




A2158

Titanium iris
repositor



A2160

Hirschman lens
spatula,
0.5mm wide



A2162

Jaffe lens spatula,
0.5mm wide



A2170

Moorfields
cystotome



A2184

Daviel expressor
spoon



Hooks, Probes and Chops

A2196

Sharp spoon
(Curette), size A



A2198

Sharp spoon
(Curette), size B



A2200

Sharp spoon
(Curette), size C



A2202

Sharp spoon
(Curette), size D



A2204

Sharp spoon
(Curette), size E



A2206

Sharp spoon
(Curette), size F



Hooks, Probes and Chops

A2208

Lang dissector and
curette



A2210

Snellen lens loop



A2211

Bell Taylor vectis



A2238

Graefe strabismus
hook, miniature



A2240

Graefe strabismus
hook, small



A2242

Graefe strabismus
hook, medium



A2244

Graefe strabismus
hook, large



Hooks, Probes and Chops

A2246

Jameson strabismus
hook



A2247

Jameson hook,
extended model



A2248

Chavasse strabismus
hook



A2250

Moorfields
strabismus hook



A2252

Arruga lens
expressor,
right



A2254

Arruga lens
expressor,
left



A2264


Guthrie iris hook,
small sharp points



Hooks, Probes and Chops

A2267

Kilner (Gillie) hook



A2268

Rollet lacrimal sac hook



A2269

Rollet lacrimal sac hook, miniature



A2277

Titanium iris hook, left



A2278

Titanium iris hook, right



A2280


Graether collar button retractor



Hooks, Probes and Chops


A2286

Hirschman iris hook



A2289

Peterborough iris
expanding hooks
(1 pair)



A2290

Kuglen IOL
hook



A2294

Fenzl IOL hook, blunt
0.2mm tip




A2297

Lester IOL hook,
0.2mm tip



A2299

Osher "Y"
hook
(horizontal tips)



A2300

Clayman IOL hook,
angled 0.2mm tip,
titanium



Hooks, Probes and Chops

A2301

Mushroom manipulator, titanium



A2302

Double ended mushroom manipulator / phaco chop, titanium




A2302E

Entabi Double ended Phaco Manipulator with mushroom




A2303

Kuglen IOL hook, titanium



A2304

Sinskey IOL hook, straight 0.2mm tip, titanium




A2305

Maltzman Fenzl hook, 0.15mm tip



A2306

Maltzman Fenzl hook, 0.15mm tip, titanium



Hooks, Probes and Chops

A2307

Hirschman spatula,
titanium




A2308

Hirschman hook,
titanium



A2310

Bechert
nucleus rotator
(horizontal tips)




A2310S

Sung manipulator
(for Classic Malyugin
Ring)





A2311

Bechert
nucleus rotator,
titanium
(horizontal tips)





A2312

Clayman nucleus
rotator





A2313

Maloney double-
ended nucleus
manipulator
(horizontal tip)



A2313A

Maloney double-
ended nucleus
manipulator
modified



Hooks, Probes and Chops

A2316

Phaco chop nucleus divider, titanium



A2316MG

Double Ended Mushroom Manipulator and Greens chop, titanium




A2316G

Green nucleus chop, titanium





A2316R

Phaco paddle, titanium




A2316RM

Phaco micro paddle, titanium




A2316RS

Nucleus splitter, titanium



A2316T

Phaco triple chop nucleus divider, titanium



A4482

Microfinger, titanium



A4482L

Microfinger, titanium for left-handed surgeons

Irrigation / Aspiration

Intertip I/A System

A3320

Intertip I/A
handpiece



A3330

Intertip,
straight



A3332

Intertip,
curved




A3334

Intertip,
45°



A3336

Intertip,
90°




A3340

Intertip, binkhorst



A3342

Intertip, straight,
with silicone sleeve



A3344

Intertip luer adapter



A3346

Replacement o-rings

Pack of 15



A3349

Intertip
flush kit



A3345

Silicone sleeves
For A3342


Pack of 5



Bimanual I/A Handpieces

A3370

Bimanual handpiece,
straight, open end,
21G irrigator




A3372

Bimanual handpiece,
curved, open end,
21G irrigator




A3374

Bimanual handpiece,
straight, closed end,
21G aspirator, 0.3mm
side port




A3376

Bimanual handpiece,
curved, closed end,
21G aspirator, 0.3mm
side port



A3378

Bimanual handpiece,
straight, closed end,
21G irrigator, 0.5mm
twin side ports



For our range of sterile Bimanual Handpieces please refer to pages 135 to 136.




Irrigation / Aspiration

Bimanual I/A Handpieces


A3380

Bimanual handpiece,
curved, closed end,
21G irrigator, 0.5mm
twin side ports



A3382

Bimanual handpiece,
curved, closed end,
23G aspirator, 0.3mm
side port




A3384

Bimanual handpiece,
curved, closed
end, 23G aspirator,
capsule polisher,
0.3mm
side port




A3386

Bimanual handpiece,
straight, closed end,
23G irrigator, 0.4mm
twin side ports



A3388


Bimanual handpiece,
curved, closed end,
23G irrigator, 0.4mm
twin side ports



Specula


A3500

Barraquer speculum,
adult size




A3500S

Barraquer sliding
speculum, adult size



A3500T

Barraquer speculum,
adult size, titanium




A3502

Barraquer speculum,
child size



A3503

Barraquer speculum,
premature size




A3505

Premature Baby
Speculum Thin
Gauge




A3513

Jaffe lid
retractor, child, 9mm
blades
(1 pair)




A3514

Jaffe lid
retractor, adult,
14mm blades
(1 pair)




A3515

Kratz Barraquer
speculum for phaco



A3515T


Kratz
Barraquer speculum
for phaco,
titanium



Specula


A3516

Kratz Barraquer
temporal incision
speculum for phaco




A3518

Kratz Barraquer
sliding speculum



A3516T

Kratz Barraquer
temporal incision
speculum for phaco,
titanium



A3520

Kratz Barraquer
adjustable speculum,
screw action



A3521

Kratz Barraquer
adjustable temporal
incision speculum,
screw action




A3523

Kratz Barraquer
adjustable speculum,
child size, screw
action



A3526

Reversible
adjustable speculum,
screw action



Specula

A3530

Clarke speculum,
adult, right



A3532

Clarke speculum,
adult, left



A3534

Clarke speculum,
child, right



A3536

Clarke speculum,
child, left



A3537

Clarke speculum,
infant, right



A3538

Clarke speculum,
infant, left



A3550

Lang speculum,
adult, right



A3552

Lang speculum,
adult, left



Specula

A3554

Lang speculum,
child, right



A3556

Lang speculum,
child, left



A3557

Lang speculum,
infant, right



A3558

Lang speculum,
infant, left



A3560

Lister speculum,
adult, right



A3562

Lister speculum,
adult, left



A3564

Lister speculum,
child, right



A3566

Lister speculum,
child, left



Specula

A3574

Cook speculum,
infant, 7mm blades



A3572

Cook speculum,
child, 11mm blades



A3576

Cook speculum,
adult, 15mm blades



A3590

Sauer speculum,
child, 10mm blades




A3592

Sauer speculum,
infant, 6mm blades



A3594

Sauer speculum,
premature, 4mm
blades



A3602

Muller retractor, 3 x 3
prong



A3624

Maumanee Park
speculum



Retractors and Rings

A3630

Desmarres retractor,
12mm blade



A3632

Desmarres retractor,
14mm blade



A3634

Desmarres retractor,
16mm blade



A3636

Desmarres retractor,
18mm blade



A3638

Desmarres retractor,
infant, 8mm blade



A3650

Fison retractor




Retractors and Rings

A3692


Flieringa ring set

Pack of 11



A3675

Flieringa ring
12mm




A3683

Flieringa ring
18mm

A3677

Flieringa ring
13mm




A3684

Flieringa ring
20mm

A3679

Flieringa ring
15mm




A3686

Flieringa ring
21mm

A3680

Flieringa ring
16mm





A3688

Flieringa ring
22mm

A3681

Flieringa ring
17mm




A3690

Flieringa ring
23mm

A3682

Flieringa ring
18mm



Needleholders

Titanium Micro Needleholders

A3850

Altomed titanium
micro needleholder,
standard, curved
115mm




A3852

Altomed titanium
micro needleholder,
standard, straight
115mm




A3860

Altomed titanium
micro needleholder,
fine, curved
115mm



A3862

Altomed titanium
micro needleholder,
fine, straight
115mm




N.B. Our titanium needleholders are also available with lock, upon request.
Suffix with 'L' eg A3850L

Needleholders

Castroviejo Micro Needleholders


A4000

Castroviejo micro
jaw needleholder,
straight
140mm




A4002

Castroviejo micro
jaw needleholder,
curved
140mm




A4004

Castroviejo micro
jaw needleholder,
straight, with lock
140mm



A4006

Castroviejo micro
jaw needleholder,
curved, with lock
140mm



Catford Micro Needleholders

A4011

Catford micro
needleholder,
straight
95mm



A4013

Catford micro
needleholder,
straight, with lock
95mm



Needleholders

Catford Micro Needleholders

A4015

Catford micro
needleholder, curved
95mm



A4017

Catford micro
needleholder,
curved, with lock
95mm



A4021

Catford micro
needleholder,
straight
128mm



A4023

Catford micro
needleholder,
straight, with lock
128mm




A4025

Catford micro
needleholder, curved
128mm



A4027

Catford micro
needleholder,
curved, with lock
128mm



Needleholders

Troutman Needleholders

A4061

Troutman
needleholder, curved
115mm




A4062

Troutman
needleholder,
curved, with lock
115mm



A4065

Troutman
needleholder,
modified with finer
extended jaws
120mm



Barraquer Needleholders

A4073

Barraquer micro
needleholder,
straight
120mm



A4075

Barraquer micro
needleholder,
straight, with lock
120mm



A4077

Barraquer micro
needleholder, curved
120mm



Needleholders

Barraquer Needleholders


A4079

Barraquer micro
needleholder,
curved, with lock
120mm



A4080

Barraquer
needleholder,
straight
135mm




A4082

Barraquer
needleholder,
straight, with lock
135mm



A4084

Barraquer
needleholder, curved
135mm



A4086

Barraquer
needleholder,
curved, with lock
135mm




Needleholders

Castroviejo Needleholders


A4112

Castroviejo
needleholder,
tungsten carbide
jaws, straight
140mm




A4113

Castroviejo
needleholder,
tungsten carbide
jaws, straight,
with lock
140mm



A4114

Castroviejo
needleholder,
tungsten carbide
jaws, curved
140mm



A4115

Castroviejo
needleholder,
tungsten carbide
jaws, curved,
with lock
140mm




Needleholders

Castroviejo Needleholders


A4161

Castroviejo
needleholder,
standard jaw,
straight
140mm




A4162

Castroviejo
needleholder,
standard jaw,
straight, with lock
140mm



A4163

Castroviejo
needleholder,
standard jaw,
curved
140mm



A4164

Castroviejo
needleholder,
standard jaw,
curved, with lock
140mm




Needleholders

Altomed Needleholder

A4088

Altomed micro
needleholder,
short jaw
128mm



Lim Needleholders


A4168CX

Lim micro
needleholder,
black finish
115mm



A4169


Lim micro
needleholder,
silver finish
115mm



Birks Needleholder

A4174

Birks micro
needleholder
115mm



Silcocks Needleholder

A4160

Silcocks
needleholder
145mm



DCR Instruments

Hajeks Punch

130mm shaft length

A4500 2mm, upward cutting

A4502 2mm, downward cutting

A4505 3mm, upward cutting

A4507 3mm, downward cutting

A4510 4mm, upward cutting

A4512 4mm, downward cutting



Images not to scale

Citelli Punch

75mm shaft length

A4528 2mm, upward cutting


A4529 2mm, downward cutting

A4530 3mm, upward cutting

A4531 3mm, downward cutting

A4532 4mm, upward cutting

A4533 4mm, downward cutting






Images not to scale

Malhotra Endonasal Punch

The new Malhotra side cutting endonasal DCR punch has been specially designed to allow easy dismantling and processing. (See instructions below)

190mm shaft length


- A4544** 45° 3mm, left cutting
- A4545** 45° 4mm, left cutting
- A4546** 45° 3mm, right cutting
- A4547** 45° 4mm, right cutting
- A4548** 90° 3mm, left cutting
- A4549** 90° 4mm, left cutting
- A4550** 90° 3mm, right cutting
- A4551** 90° 4mm, right cutting




Images not to scale

To Dismantle:

1. Close the handle and then pull back the release trigger.





2. Carefully slide back the moveable arm and separate from handle.



Important! Ensure the two parts of the device are kept together and are not separated from one another. Each arm is custom made for each handle.

To Reassemble:

After cleaning, lubricating and sterilising, reassemble the punch by aligning the Slip Plate with the Slip Groove and insert the Retaining Pin into the Retaining Slot. Carefully close the handle and push back the Release Trigger into the closed position.



DCR Instruments

Sphenoidal Punch

A4536

Sphenoidal punch,
compound action
200mm



A4537

Malhotra Endonasal
DCR nibbler delicate
200mm



Images not to scale

DCR Trephines

A4538

Arruga Rycroft DCR
trephine, coarse
teeth,
with pin
190mm



A4539

Arruga Rycroft DCR
trephine,
fine teeth,
without pin
190mm



Images not to scale

A4570

Bone chisel,
3mm tip

**A4572**

Bone chisel,
4mm tip

**A4574**

Bone chisel,
5mm tip

**A4580**

Bone gouge,
3mm tip

**A4582**

Bone gouge,
4mm tip

**A4584**

Bone gouge,
5mm tip

**A4586**

Mallet
210mm



Images not to scale

A4595

Lang sharp spoon,
double ended



DCR Instruments

A4599

Freer
elevator
180mm



Images not to scale

A4600

Traquair periosteal
elevator
160mm



Images not to scale

A4601

Traquair lacrimal
rasp



A4602

Stallards dissector,
blunt




A4604

Rollet
rugine



A4605


Howarth
elevator
210mm



Images not to scale

A4606

Rollet rugine / elevator
double ended
160mm



Images not to scale

The Crawford Lacrimal Intubation System

This is a simple procedure designed by the late J.S. Crawford, M.D, which speeds up the process of inserting a lacrimal stent with less trauma and tearing of the nasal floor.

Indications for lacrimal intubation:

1. Where probing of the lacrimal system has not cured the tearing.
2. Where there are definite obstructions of the lacrimal system that would close if they were not kept open.
3. Repair of injuries involving cut canaliculi.
4. During a dacryocystorhinostomy.

28-0184 Crawford intubation system with suture

Two flexible stainless steel probes (0.40mm diameter wires) with an olive shaped tip at one end, with the other ends joined by a length of silicone tubing (outside diameter 0.64mm, inside diameter 0.30mm). This version of the Crawford lacrimal intubation set contains a 6-0 braided silk suture within the silicone tube lumen.

28-0185 As above but without suture

28-0174 Crawford II intubation system with suture

Two flexible stainless steel probes (0.40mm diameter wires) with an olive shaped tip at one end, with the other ends joined by a length of silicone tubing (outside diameter 0.94mm, inside diameter 0.51mm). This version of the Crawford lacrimal intubation set contains a 6-0 braided silk suture within the silicone tube lumen.

28-0175 As above but without suture

28-0186 Crawford retrieval hook

Used to pick up the probe (after insertion) under the inferior turbinate in the nose.

28-0189 Anderson-Hwang grooved director


This grooved director greatly simplifies retrieval of the olive tip probe. The grooved director has the following advantages:

Minimizes trauma and bleeding of the nasal mucosa.
Greatly simplifies the intubation technique by trapping the olive tip probe.
Allows infraction of the inferior turbinate if required.

The grooved director also simplifies placement of tubes during a DCR. The olive tip of the Crawford set may be engaged either under direct visualization or with an endoscope.

28-0190 Crawford tubing stripper

Used to cut and remove excess tubing from the silk suture.




Note: Intubation sets are supplied sterile and are Single Use products. The hook, stripper and grooved director are intended to be sterilized at the hospital and can be reused provided they are not damaged.

CDCR Tubes

Our range of implantable pyrex glass tubes is intended for use in CDCR procedures. There is a wide variety of styles, with a selection of flange diameters and tube lengths, each tube being manufactured to our customers specifications.
Please refer to the Altomed price list for guidance.


Lester Jones Tubes, Straight and Angled



Manufactured from 2.2mm diameter pyrex borosilicate glass tubing with a 3mm, 3.5mm or 4mm flange.
Manufactured with a bevel at the nasal end.

Order code LJ_____

Gladstone Putterman Tube



Similar to a Lester Jones tube but has a second cone-shaped flange 4-6mm behind outer 3mm, 3.5mm or 4mm flange. This has been developed to prevent tube migration. When ordering quote distance to 2nd flange.

Order code GP_____

Callahan (Cox) Tube



Manufactured from 3mm pyrex tubing, normal flange diameters 4.5mm or 5.5mm, has no nasal end bevel.

Order code CC_____

Double Collar Tube



This is a straight Lester Jones tube with a flange at both ends, supplied with either 3mm, 3.5mm or 4mm flange and minimum length of 12mm.

Order code DC_____

Ordering Guide


LJ3509 Lester Jones tube, 3.5mm flange, 9mm length.

LJ4015AHF Lester Jones tube, 4mm flange, 15mm length, angled, with suture hole, frosted.

DC3017F Double collar tube, 3mm flange, 17mm length, frosted.

- A11908S** Polyethylene tube, sterile (pack of 5).
- A11930** Plastic storage tray with labelled compartments (not autoclavable).
- A11932** Blue anodized aluminium autoclavable sterilisation tray.
- A11934** Blue anodized aluminium identification scale for tubes.
- A11936** Pack of 5 Jones nylon cleaning rods (not autoclavable).
- A11937** Lester Jones gold dilator No.I 1.7mm and 1.9mm arms.
- A11937A** Lester Jones gold dilator No.II 2.2mm and 2.4mm arms.
- A11937B** Lester Jones gold dilator No.III 2.7mm and 2.9mm arms.

Polyethylene Tube



Jones Gold Dilators



In response to the mechanical difficulties observed during conjunctivo-dacryocystorhinostomy, a set of instruments was developed by C.W. Cox Jr, to make the procedure technically easier with less trauma to the patient, and successful functioning of the tube more likely. With the cooperation of oculoplastic surgeons, notably Dr Alston Callahan, the modern Cox System has been developed.

COX2

Measuring gauge

**COX03**

Trocar

**COX04**

Teflon pusher

**COX05**

Insertion rod, small

**COX06**

Insertion rod, large

**COX07**

Flat ended rod for measuring tube lengths

**A11940**

Parkin tube introducer



Full instructions regarding the use of the Cox System are available on request, or can be downloaded from our website.

Depressors, Punches and Calipers

A4800CX

Schepens localiser, size 8
thimble (small)



A4802CX

Schepens localiser, size 10
thimble (medium)



A4804CX

Schepens localiser, size 12
thimble (large)



A4810CX

Gass scleral marker, size 12
thimble, 1.5mm tip



A4813

Packo Combined
Depressor and
Marker



A4814

Schocket scleral
depressor




A4815

Neonate scleral
depressor (Schocket
type for infants)




A4818

Keratech astigmatic
rule (reusable and
autoclavable)



A4819


Maloney
keratoscope



Depressors, Punches and Calipers

A4830

Kelly punch,
1.0mm



A4842

Membrane punch,
0.7mm



A4900

Castroviejo caliper, straight,
15mm spread




A4902

Castroviejo caliper, curved,
20mm spread



A4904

Rule,
stainless steel
150mm



A4906

Giessen rule,
30mm



A4908

Sclerotomy caliper,
titanium
(3.0, 3.5, 4.0mm)



Images not to scale

Pressure Reducers and Corneal Burrs

A4949

Honan balloon pressure reducer with relief valve

A4950

Spare headband for Honan balloon

A4951D

Disposable bellows and headband for Honan balloon (box of 12)



A4955

Algerbrush corneal rust ring remover



A4959

Single Use burr for Algerbrush (box of 5)



A4960


Single Use burr for Algerbrush (box of 25)



Sterile

PD-25RFD


AlgerBrush 2.5mm Pterygium Diamond Burr Unsterile (each)



Bipolar


A5002

McPherson
bipolar forceps,
straight,
2 pin fitting
90mm




A5004

McPherson
bipolar forceps,
angled,
2 pin fitting
90mm




A5006

Jewellers
bipolar forceps,
2 pin fitting
100mm




A5040

Bipolar forceps,
straight,
BS fitting
80mm




A5042

Bipolar forceps,
straight,
BS fitting
110mm



A5044


Bipolar forceps,
angled,
BS fitting
110mm



Bipolar Forceps and Pencils


A5040P

Bipolar forceps,
straight,
2 pin fitting
80mm




A5042P

Bipolar forceps,
straight,
2 pin fitting
110mm




A5044P

Bipolar forceps,
angled,
2 pin fitting
110mm



A5050

Bipolar pencil,
18G,
straight,
2 pin fitting reusable





A5052

Bipolar pencil,
18G,
angled,
2 pin fitting reusable




A5056

Bipolar pencil,
20G,
straight,
2 pin fitting reusable



A5058

Bipolar pencil,
20/23G, tapered,
2 pin fitting reusable



For our range of sterile Bipolar instruments please refer to page 102.



Bipolar Cables

Altomed offers a wide range of bipolar cables available for electrosurgery.

Choose an Instrument fitting from the tables to the right followed by a supporting machine fitting, to specify the product code.

* 2 Plug 22mm block compatible with Infiniti machine.

† 2 Plug 28.6mm block compatible with Accurus, Constellation, Centurion, Aesculap, Valleylab machines.

Instrument Fittings	Machine Fittings	Code
 2 pin	2 plug 22mm block	A5080-22
	2 plug 28.6mm block	A5080-28.6
	Half round	A5081
	British Standard / Millennium	A5098
	Oertli Blue Ring	2058-72
	Oertli Red Ring	2063-72

Instrument Fittings	Machine Fittings	Code
 British Standard / Millennium	British Standard / Millennium	A5084
	Half round	A5085
	2 plug 22mm block	A5086-22
	2 plug 28.6mm block	A5086-28.6
	Oertli Blue Ring	2058-71
	Oertli Red Ring	2063-71

Machine Fittings



British Standard / Millennium



2 plug 22mm block *



Oertli Blue Ring



Half Round



2 plug 28.6mm block †



Oertli Red Ring

N.B. Other configurations of Bipolar Cable can be supplied upon request.

Forceps

A5100

Halstead Mosquito
artery forceps,
straight
120mm



A5102

Halstead Mosquito
artery forceps,
curved
120mm



A5110

Halstead Mosquito
artery forceps,
straight
100mm



A5112

Halstead Mosquito
artery forceps,
curved
100mm



Forceps

A5114

Spencer Wells artery
forceps
180mm



A5116

Ramsey sponge
holder
180mm



A5117

Tilley nasal forceps
150mm



Images not to scale

Forceps

A5119

Cross action towel clips
90mm



A5120

Backhaus towel clips
90mm



A5124

Ball & socket towel clips
90mm



A5130

Bulldog clip, straight
25mm



A5132

Bulldog clip, straight
40mm



A5140

Bulldog clip, curved
25mm



A5142



Bulldog clip, curved
40mm



Forceps


A3800

Straight titanium
micro forceps
90mm




A3801

Straight titanium
long handled micro
forceps
120mm




A3802

Straight titanium
micro forceps,
0.12mm 1/2 teeth
90mm




A3804

Straight titanium
micro forceps,
0.25mm notched
90mm





A3810

Angled titanium
micro forceps
90mm



A3812

Angled titanium
micro forceps,
0.12mm 1/2 teeth
90mm



A3814

Angled titanium
micro forceps,
0.25mm notched
90mm



Forceps

A3820

Colibri titanium
micro forceps
90mm



A3822

Colibri titanium
micro forceps, 0.12
mm
1/2 teeth
90mm





A3824

Colibri titanium
micro forceps,
0.25mm notched
90mm



A3830

Curved titanium
micro forceps
90mm





A3831

Curved titanium
long handled micro
forceps
120mm




A3832

Curved titanium
micro forceps,
0.12mm
1/2 teeth
90mm



A3834


Curved titanium
micro forceps,
0.25mm notched
90mm



Forceps


A3840

Swan necked
titanium micro
forceps
90mm




A3842

Swan necked
titanium micro
forceps, 0.12mm
1/2 teeth
90mm



A3844


Swan necked
titanium micro
forceps,
0.25mm notched
90mm



Forceps

A5400

Arruga intracapsular
forceps
90mm




A5416

Castroviejo capsular
forceps
100mm




A5450

Barraquer cilia
forceps
105mm



A5452

Bennett cilia forceps
70mm



A5454

Beer cilia forceps
90mm



A5462


Whitfield cilia
forceps
90mm



Forceps


A5498

Schulenburg
Western forceps
110mm




A5500

Moorfields forceps
110mm




A5500F

Moorfields forceps
fine tips
106mm




A5502

Listers conjunctival
forceps
100mm




A5506

Graefe dissecting
forceps,
serrated tips
100mm




A5510

Adsons dissecting
forceps,
serrated tips
125mm



A5512


Adsons dissecting
forceps,
1/2 teeth
125mm



Forceps

A5542

Castroviejo tying
forceps,
0.12mm teeth
100mm




A5558

St. Martins suture
forceps, micro,
0.12mm teeth
85mm




A5559

St. Martins suture
forceps,
0.3mm teeth
85mm




A5560

Jayles tying forceps,
0.3mm teeth
115mm



A5562


Jayles tying forceps,
micro,
0.12mm teeth
100mm



Forceps


A5580

Hoskins forceps,
No.18,
0.5mm tips,
Silver titanium




A5582

Hoskins forceps,
No.19,
0.3mm tips,
Silver titanium





A5586

Hoskins forceps,
No.22,
0.3mm tips,
Silver titanium



A5588

Hoskins forceps,
No.28,
0.25mm tips,
Silver titanium





A5590

Hoskins forceps,
No.30,
0.25mm tips,
Silver titanium



A5592


Hoskins forceps,
No.32,
0.3mm tips,
Silver titanium



Forceps

A5656

Barraquer colibri
corneal tying
forceps,
0.12mm teeth at 90°
80mm




A5657CX

Barraquer colibri
corneal tying
forceps,
0.12mm teeth at 90°,
blackened
80mm




A5660

Barraquer colibri
corneal tying
forceps,
0.12mm teeth at 45°
80mm




A5663CX

Barraquer colibri
corneal tying
forceps,
0.12mm teeth at 45°,
blackened
80mm



A5674

Polack double
corneal forceps,
2 x 0.12mm teeth
75mm



A5674F



Polack double
corneal forceps,
extra fine teeth
75mm



Forceps


A5682

McPherson corneal forceps,
0.12mm 1/2 teeth,
straight
85mm




A5686

Kelman-McPherson corneal forceps,
0.12mm 1/2 teeth,
angled
85mm




A5688

Bonn corneal forceps,
0.12mm 1/2 teeth,
straight
85mm




A5804

Thaller conjunctival forceps
90mm




A5806

Elschnig fixation forceps
0.12mm
1/2 teeth
110mm



A5812

Dastoor forceps
110mm



A5828


Barraquer fixation forceps,
multiple teeth
100mm



Forceps


A5832

Francis fixation
forceps,
2/3 teeth
110mm




A5833

Francis fixation
forceps,
3/4 teeth
110mm




A5834

French fixation
forceps,
2/3 teeth
110mm




A5835

French fixation
forceps,
3/4 teeth
110mm




A5854

Graefe iris forceps,
1/2 teeth
75mm



A5855

Barraquer iris
forceps,
serrated tips
70mm



A5857


Barraquer iris
forceps,
1/2 teeth
70mm



Forceps


A5884

Silcocks iris forceps,
1/2 teeth
100mm




A5886

Silcocks iris forceps,
serrated tips
100mm




A5900

Jewellers forceps,
No.1
110mm



A5908

Jewellers forceps,
No.5
110mm




A5916

Jewellers forceps,
No.7
110mm





A5926

Troutman superior
rectus forceps,
1/2 teeth
110mm



A5927


Moody fixation
forceps, with lock,
1/2 teeth
110mm



Forceps

A5937

IOL folding forceps,
silver titanium
110mm





A5938

IOL insertion forceps,
silver titanium
110mm




A5940

Soft lens folding
and holding forceps,
titanium
110mm



A5944

Soft lens insertion
forceps,
titanium
100mm





A5944L

Soft lens insertion
forceps,
with lock,
titanium
100mm



A5945

Modified soft lens
folding and holding
forceps,
titanium
110mm



Forceps

A5958

Clayman IOL forceps,
right hand
120mm




A5960

Clayman IOL forceps,
left hand
120mm




A5966

Shepard IOL forceps
120mm




A5967

T-end IOL forceps,
titanium
115mm




A5969

Clayman IOL forceps,
titanium
115mm



A5970

IOL Cutter
130mm




Images not to scale

Forceps


A5976

McPherson tying
forceps,
straight,
5mm platform
85mm




A5978

McPherson tying
forceps,
angled,
5mm platform
85mm



A5980

Kelman McPherson
forceps,
angled, 8mm jaw,
5mm platform
82mm



A5981

Kelman McPherson
forceps,
angled,
8mm platform
100mm




A5983CX

Round body micro
forceps,
straight, blackened
100mm



A5985CX

Round body micro
forceps, straight,
0.12mm
1/2 teeth,
blackened
100mm



Forceps

A5987CX

Round body micro
forceps,
curved,
blackened
100mm





A5989CX

Round body micro
forceps,
curved,
0.12mm teeth,
blackened
100mm




A6011

Mermound Deep
Sclerotomy forceps




A6014

Sheets McPherson
forceps,
10mm jaw,
serrated tips
100mm



A6017F

Capsulorhexis
forceps,
X action
11mm jaw,
Serrated,
titanium
100mm



A6020

Utrata Capsulorhexis
forceps,
11mm jaw,
silver titanium
100mm



A6021

Utrata Capsulorhexis
forceps,
11mm curved jaw,
silver titanium
100mm



N.B. Capsulorhexis tip details shown under high magnification

Forceps

A6022

Utrata Capsulorhexis
forceps, round
handle,
titanium
100mm




A6023

Utrata Capsulorhexis
forceps,
round handle,
curved jaw,
titanium
100mm




A6024

Utrata Capsulorhexis
forceps,
titanium
85mm



A6025

Utrata Capsulorhexis
forceps,
curved jaw,
titanium
85mm




A6027

Capsulorhexis
forceps, extra fine,
curved jaw round
handle
100mm



A6028

Coventry Posterior
Capsulorhexis
forceps,
titanium



A6029

Capsulorhexis
forceps,
micro serrated tips
85mm




N.B. Capsulorhexis tip details shown under high magnification.

Forceps

A6030


Capsulorhexis
forceps,
21G
140mm



Images not to scale

A6034


Capsulorhexis
forceps,
23G
140mm



Images not to scale


A6162

Rycroft suture
forceps
100mm




A6169

Max fine tying
forceps,
curved
100mm





A6170

Birks (Tennant)
micro tying forceps,
straight
100mm




A6170T

Birks (Tennant)
micro tying forceps,
straight,
titanium
100mm



A6172

Birks (Tennant) micro
tying forceps,
curved
100mm



A6172T

Birks (Tennant)
micro tying forceps,
curved,
titanium
100mm



Forceps

A6100

Jameson strabismus
forceps,
seven prong, left
100mm



A6102

Jameson strabismus
forceps,
seven prong, right
100mm



A6108

Prince strabismus
forceps,
three prong, left
100mm



A6110

Prince strabismus
forceps,
three prong, right
100mm




A6384

Watzke forceps for
silicone sleeve
120mm



A6388

Watzke forceps for
silicone sleeve, with
fixing screw
120mm



Lid Surgery

A6300

Berke ptosis forceps
110mm



A6304

Erhardt lid
forceps
90mm



A6312

Stallard lid spatula
(pair of large and
small)



A6316

Macallans lid spatula
90mm



A6370

Wrights fascia needle
130mm



A6370P

Wrights fascia
needle,
child size
130mm



Images not to scale

Lid Surgery

A6320

Desmarres Entropian
forceps, left



A6322

Desmarres Entropian
forceps, right




A6330

Snellen Entropian
forceps, left




A6332

Snellen Entropian
forceps, right



A6342

Wilde Entropian
forceps



Lid Surgery

A6344

Desmarres chalazion
forceps,
16mm



A6345

Desmarres chalazion
forceps,
20mm



A6346

Desmarres chalazion
forceps,
26mm



A6347

Desmarres chalazion
forceps,
31mm



A6350

Ayers tarsal cyst
forceps



A6352


Cross action tarsal
cyst forceps



Lid Surgery

A6354

Birmingham tarsal cyst forceps



A6356

Greene tarsal cyst forceps



A6358

Lambert tarsal cyst forceps



A6376

Protective eye shell, autoclavable (includes introducer)



Corneal Markers and Global Fixation

A6600

Thornton
fixation ring,
16mm diameter



A6602

Thornton
fixation ring,
with gap,
16mm diameter



A6604

Thornton
fixation ring,
with gap
and swivel,
16mm diameter



A6605

Thornton
fixation ring,
with gap
and swivel,
13mm diameter



A6640

Mendez degree
gauge




A6642

Bores axis marker



A6644


Thornton limbal
incision rule



Corneal Markers and Global Fixation


A6692

Osher Neuman
corneal marker,
eight lines




A6693

Osher Neuman
corneal marker,
sixteen lines



A6694

Osher Neuman
corneal marker,
twelve lines



Corneal Trephines and Punches

A6910

Elliot trephine handle



A6912

Elliot trephine blade,
1.5mm
set of two



A6914

Elliot trephine blade,
2.0mm
set of two



For our range of sterile Trephine blades and cutting blocks please refer to page 99.



Enucleation and Orbit

A7060

Trelat spatula
150mm




A7062

Elliot spatula
150mm




A7064

Mule evisceration
spatula
115mm



A7074

Wright orbital
forceps
120mm



Vitreoretinal Instruments

A7662

Scleral plug forceps



For disposable Scleral plugs please refer to page 107.



Volk Lenses

Classic Indirect Bio Lenses

A7804

20D Lens



A7804

V20LCD10

20D Lens Sterile, Single Use
Box of 10



V20LCD10



Box of 10

A7819

28D Lens



A7819

V28LCD10

28D Lens Sterile, Single Use
Box of 10



V28LCD10



A7812

30D Lens



A7812



A7826

A7826

Pan Retinal 2.2 Lens

Images not to scale

Digital Indirect Bio Lenses

VDGTLCM

Digital Clearmag Lens



VDGTLCM

VDGTLCF

Digital Clearfield Lens



VDGTLCF

Images not to scale

Slit Lamp Lenses

A7799

Super 66 Lens



A7799



A7828

A7828

Superfield NC Lens

A7824

78D Lens



A7824



A7822

A7822

90D Lens



A7792

A7792

Super Vitreofundus Lens

Digital Slit Lamp Lenses

VDGTL1

Digital 1.0x Imaging Lens



VDGTL1



VDGTLHM

VDGTLHM

Digital High Mag Lens

VDGTLWF

Digital Wide Field Lens



VDGTLWF

Images not to scale

YAG Laser Lenses

A7827C
Capsulotomy Laser Lens



A7827C



VCAPSD10

VCAPSD10
Capsulotomy Laser Lens
Sterile Single Use
Box of 10

A7827R
Iridectomy Laser Lens



A7827R



VIRIDD10

VIRIDD10
Iridectomy Laser Lens
Sterile Single Use
Box of 10

A7827M
Magplus Iridectomy Lens



A7827M

Images not to scale

Indirect Contact Laser Lenses

A7839
Superquad 160 Lens



A7839



A7840

A7840
Quadraspheric Lens

A7841
Area Centralis Lens



A7841



A7842

A7842
Trans Equator Lens

VHRC
High Resolution Centralis Lens



VHRC

VHRWF
High Resolution Wide Field Lens



VHRWF

Images not to scale

Volk Lenses

Gonio Lenses

VG1

One Mirror Gonio Lens *

VG3

Three Mirror Gonio Lens *

V3MIRD10

Three Mirror Gonio Lens
Sterile Single Use
Box of 10

VG4HAN2

Four Mirror Gonio Lens, on two
position handle

VG4

Four Mirror Gonio Lens *

V4MIRD10

Four Mirror Gonio Lens
Sterile Single Use
Box of 10

VSLTD10

SLT Lens
Sterile Single Use
Box of 10

* 'No-flange' versions are available
by suffixing NF



VG1



VG3



VG4



VG4HAN2



V3MIRD10



V4MIRD10



VSLTD10



Images not to scale

Accessories

A7846

Precision Optical Lens Cleaner

Box of 24



New Volk
Single-Lens Cases
features a sleek
and modern
functional design.

Microwash Trays and Pillows

Microwash Sterilization Trays

A9304

Microwash tray
190 x 75 x 40mm
(1/8 DIN)

A9306

Microwash tray
230 x 125 x 40mm
(1/4 DIN)

A9310

Microwash tray
230 x 255 x 40mm
(1/2 DIN)

A9311

Microwash tray
340 x 255 x 40mm
(3/4 DIN)

A9312


Microwash tray
445 x 255 x 40mm
(1 DIN)

A9313


Microwash tray
340 x 250 x 80mm
(3/4 DIN
double deck)

A9314

Microwash tray
445 x 255 x 80mm
(1 DIN, double deck)



Microwash trays enable microsurgical instruments to be washed and sterilized using modern automated methods. Microwash trays accommodate instruments securely, preventing damage during handling, whilst also allowing fluid spray penetration. The trays range from 75mm up to 460mm, with the option of single or double layers. The internal components are custom made to secure each individual instrument.



Microwash trays are supplied with laminated placement sheets which enable easy identification of instruments.

Operating Table Pillows

A7880

Rubens pillow,
adult



Sterilization Trays

Polymer Sterilization Trays with Silicone Rubber Finger Mats, Latex free

A9201

150mm x 50mm x 19mm

A9202

150mm x 50mm x 19mm
Silicone bars

A9203

150mm x 75mm x 19mm

A9210

280mm x 180mm x 18mm

A9220

380mm x 255mm x 18mm

A9222

380mm x 255mm x 38mm
with insert, space below

A9230

450mm x 250mm x 18mm

A9232

450mm x 250mm x 38mm
with insert, space below

A9239

254mm x 152mm x 18mm

A9240

254mm x 152mm x 38mm

A9241

254mm x 152mm x 38mm
with insert, space below

A9242

254mm x 152mm x 38mm
with insert, mat below


A9246

380mm x 255mm
with insert, space below

A9248

380mm x 255mm
with insert, mat below

All dimensions are internal



Surgeons Stools and Q-Rinse

Operating Stools, Hydraulic

A11020

Waterfall seat
(pictured, left)

A11022

Saddle seat

A11024

Waterfall seat, small
armrests

A11026

Saddle seat, small armrests

A11028

Waterfall seat, large
armrests

A11029

Saddle seat, large
armrests
(pictured, right)



Q-Rinse System

A11070

Q-Rinse Automated Rinse
System

Includes:

- 1x Q-Rinse machine
- 1x Wire brush
- 1x Power adapter and lead
- 1x Autoclavable bottle
assembly and tubing
- 8x Adapters
- 1x Instruction book

A11072

Q-Rinse Bottle Assembly –
Complete

A11076

Q-Rinse Wire Brush Single
Use Non-Sterile



Bottle Parts

A11044A

Q-Rinse Bottle Only
No Cap or Tubing

A11044C

Q-Rinse Bottle Cap Gaskets
Pack of 4

Further accessories and spare parts are available to order.
Please contact customer service for details.



We repair all makes

Altomed is able to repair any brand of ophthalmic surgical instrument or diamond knife.

We can also repair almost any type of phacoemulsification handpiece.

Once refurbished, your old handpieces will be as good as new and the repair will be covered by our own warranty.

Free estimates are available. Ask for 'Repairs' on 0191 5190111 or email: repairs@altomed.com
Helping you get the most from valuable instruments.



Sterile Single Use

Corneal Trephines and Punches

Hessburg-Barron Vacuum Trephines

A7030

6.50mm

A7032

7.00mm

A7033

7.25mm

A7034

7.50mm

A7035

7.75mm

A7036

8.00mm

A7037

8.25mm

A7038

8.50mm

A7039

8.75mm

A7040

9.00mm



The Hessburg-Barron trephine can provide either a complete or partial keratotomy for keratoplasty procedures. It adheres to the cornea during the process through the creation of a vacuum, holding it in its natural state. This allows for a perfect, perpendicular cut of the recipient well. Each full rotation advances the blade by 250 microns.

Barron Marking Punches

A7042

6.50mm

A7043

6.75mm

A7044

7.00mm

A7045

7.25mm

A7046

7.50mm

A7047

7.75mm

A7048

8.00mm

A7049

8.25mm

A7050

8.50mm

A7051


8.75mm

A7052

9.00mm

A7053

9.50mm



The Barron marking corneal donor punch system is designed to make an exact cut of the donor tissue for corneal transplants. The ultra-sharp cutting blade is guided by four posts to assure precise alignment. After placing the tissue in the punch, an ink pen is used to mark the donor buttons in four symmetrical locations, aiding the placement of primary sutures.

Corneal Trephines and Punches

Blade Styles

Short Trephine



Long Trephine



Short Eyetec Trephine



Long Eyetec Trephine





Size / Part No. Guide

Blade Diameter	Short Trephines	Long Trephines	Short Eyetec Trephines	Long Eyetec Trephines
2.50mm		A6816L		
3.00mm		A6818L		
3.50mm		A6820L		
4.00mm		A6822L		
4.50mm		A6824L		
5.00mm	A6826	A6826L		
5.50mm	A6828	A6828L		
6.00mm	A6832	A6832L	51-900-6.00	
6.25mm	A6834	A6834L	51-900-6.25	51-903-6.25
6.50mm	A6836	A6836L	51-900-6.50	51-903-6.50
6.75mm	A6838	A6838L	51-900-6.75	51-903-6.75
7.00mm	A6840	A6840L	51-900-7.00	51-903-7.00
7.25mm	A6842	A6842L	51-900-7.25	51-903-7.25
7.50mm	A6844	A6844L	51-900-7.50	51-903-7.50
7.75mm	A6846	A6846L	51-900-7.75	51-903-7.75
8.00mm	A6848	A6848L	51-900-8.00	51-903-8.00
8.25mm	A6850	A6850L	51-900-8.25	51-903-8.25
8.50mm	A6852	A6852L	51-900-8.50	51-903-8.50
8.75mm	A6854	A6854L	51-900-8.75	51-903-8.75
9.00mm	A6856	A6856L	51-900-9.00	51-903-9.00
9.25mm	A6858	A6858L	51-900-9.25	51-903-9.25
9.50mm	A6860	A6860L	51-900-9.50	51-903-9.50
9.75mm	A6862	A6862L		
10.00mm	A6868	A6868L		
10.50mm	A6870	A6870L		
11.00mm	A6872	A6872L		
11.50mm	A6874	A6874L		
12.00mm	A6876	A6876L		
12.50mm	A6878	A6878L		
13.00mm	A6880	A6880L		
13.50mm	A6882	A6882L		
14.00mm	A6884	A6884L		


A7012S

Teflon block for
corneal punches
A7010 and A7011
Box of 5
Sterile




A7022S

Macrae Teflon block
for Corneal punches
A7010 and A7011
Box of 5
Sterile



A7023S

Large Teflon block,
Flat base
Box of 5
Sterile




Single Use

DCR

O'Donoghue DCR Bodkins

M9290


O'Donoghue double ended DCR bodkin, 45mm, straight



Box of 5

M9291

O'Donoghue double ended DCR bodkin, 45mm, angled



Box of 5

The Crawford Lacrimal Intubation System

28-0184

Crawford intubation system with suture
Two flexible stainless steel probes (0.40mm diameter wires) with an olive shaped tip at one end, with the other ends joined by a length of silicone tubing (outside diameter 0.64mm, inside diameter 0.30mm). This version of the Crawford lacrimal intubation set contains a 6-0 braided silk suture within the silicone tube lumen.

28-0185

As above but without suture

28-0174

Crawford II intubation system with suture
Two flexible stainless steel probes (0.40mm diameter wires) with an olive shaped tip at one end, with the other ends joined by a length of silicone tubing (outside diameter 0.94mm, inside diameter 0.51mm). This version of the Crawford lacrimal intubation set contains a 6-0 braided silk suture within the silicone tube lumen.

28-0175

As above but without suture

28-0186

Crawford retrieval hook
Used to pick up the probe (after insertion) under the inferior turbinate in the nose. (reusable)



Oculoplastic Oncology

A3710

Symblepheron ring, small I.D.
12mm, O.D 21mm



A3712

Symblepheron ring, medium
I.D. 12.5mm, O.D 21mm



A3714

Symblepheron ring, large
I.D. 13mm, O.D 22mm



A7077

Damato Ruthenium Plaque
Template Small



A7078

Damato Ruthenium Plaque
Template Medium



A7079

Damato Ruthenium Plaque
Template Large



A7080

Malhotra Platinum Segment 0.2g



A7082

Malhotra Platinum Segment 0.4g



Actual size

Silicone Spheres

A7102 12mm

A7104 14mm


A7106 16mm

A7108 18mm

A7109 19mm

A7110 20mm

A7112 22mm



A7198S

Tantalum markers,
2.5mm diameter
Set of 4



Lucite Conformers

A7122 Small

A7124 Medium

A7126 Large



Silicone Conformers

A7130 Small

A7132 Medium

A7134 Large




Bipolar

Images not to scale

A5051

Bipolar pencil, 18G,
straight
120mm
Box of 10



A5053

Bipolar pencil, 18G,
angled
120mm
Box of 10




A5054

Bipolar pencil, 27G,
straight
135mm
Box of 5




A5055

Bipolar pencil, 23G,
straight
135mm
Box of 5




A5057

Bipolar pencil, 20G,
straight
135mm
Box of 10




A5063

Bipolar pencil, 25G,
straight
135mm
Box of 5




A5059

Bipolar pencil,
tapered 20G to 23G
135mm
Box of 10




A5064

Bipolar forceps,
jewellers
110mm
Box of 10



A5060 - 22

Bipolar cables,
2pin - 2 plug
22mm block
Infiniti fitting



Box of 10



Box of 10

A5060 - 28.6

Bipolar cables,
2pin - 2 plug
28.6mm block
Accurus,
Constellation or
Centurion fitting

For our range of reusable Bipolar instruments please refer to pages 57 to 59




All retinal implants are supplied in boxes of 5



Round Scleral Sponge


A7298

Style 8215




A7299

Style 502




A7289

Style 5025




A7300

Style 503




A7302

Style 504




A7304

Style 505



A7306


Style 8270



Oval Scleral Sponge


A7308

Style 501




A7310

Style 507




A7312

Style 506



A7314


Style 907



Oblong Scleral Sponge


A7286

Style 509



A7288


Style 508



Scholda Sponge


A7307

Style 8550



A7309


Style 8570



Tunnel Scleral Sponge

A7305


Style 505T




Grooved Scleral Sponge

A7250


Style 516G

**A7251**


Style 510G

**A7252**


Style 506G

**A7253**


Style 511G

**A7254**


Style 519G

**A7256**


Style 509G

**A7258**


Style 508G

**A7260**

Style 505G

**A7262**


Style 507G




Part Thickness Scleral Sponge

A7265


Style 8425

**A7266**


Style 517

**A7268**


Style 515

**A7270**

Style 510

**A7272**


Style 511




Scleral Sponge Accessories

A7280

Style 513

**A7282**


Style 514




Silicone Strips

A7398


Style 41

**A7400**

Style 42

**A7401**


Style 4050



Circling Bands


A7340

Style 40



A7370


Style 240



Silicone Lace


A7377

Style 40



A7379


Style 240



Grooved Strips


A7342

Style 20




A7344

Style 31



A7346

Style 32



A7372

Style 219




A7374

Style 220



A7376


Style 225



Sleeves


A7414

Style 270




A7416

Style 70




A7417

Style 72




A7418

Style 3084



A7419

Style 3083



A7420

Wong sleeve style 3093



A7422


Wong sleeve style 3094



Tantalum Clips


A7358

Style 50



A7392

Style 250



Tyre Implants

A7347

Style 286

**A7384**

Style 278

A7348


Style 78

**A7386**

Style 279

A7378


Style 277

**A7387**

Style 288

A7380


Style 287

**A7388**

Style 280

A7380L


Style 287WG

**A7388L**

Style 280LG

A7381


Style 275

**A7390**

Style 281

A7382

Style 276



**A7391**

Style 289

Miscellaneous Components

A7360



Style 60 Boat

**A7404**

Style 103 Meridional

A7394


Style 260 Boat

**A7406**

Style 106 Meridional

A7395


Style 135 Wedge

**A7410**

Style 112 Meridional

A7396


Style 137 Wedge

**A7412**

Style 80 Pad


A7402

Style 22 Button



A7667

Disposable
scleral plugs
(2x20G, 2x19G)



Set of 4

Flute Handles

A7669


Flute handle, passive
aspiration



Each

A7669A

Flute handle,
active or passive
aspiration



A7671/10S

Replacement
backflush tubes



Pack of 10

For our full range of Single Use cannula and handles please refer to pages 138 - 142

Infusion Cannula

A7735

20G, 2.5mm



A7736

20G, 4mm




Box of 10

Single Use

Vitreoretinal

PD.100.00

Endo illuminator,
straight, 20G
(Occutome fitting)



Box of 5

A9520A

Endo illuminator,
angled, 20G
(Occutome fitting)



Box of 10


A9520AC

Endo illuminator,
angled, 20G
(Constellation fitting)

Classic Indirect Bio Lenses

V20LCD10

20D Lens
Box of 10



V20LCD10

V28LCD10

28D Lens
Box of 10



V28LCD10

Box of 10



Images not to scale


Vitrectomy Lenses

VFD

Flat lens
36° static F.O.V.
1.02x Mag.



Plano anterior surface affords a 36° field of view of the central posterior pole and vitreous in phakic and pseudophakic eyes. This lens is ideal for photography.




Box of 10

VWFD

Widefield lens
48° static F.O.V.
0.49x Mag.



Concave anterior surface facilitates a 48° field of view when visualizing the central posterior pole and central vitreous in phakic and pseudophakic eyes.




VBCD

Biconcave lens
25° static F.O.V.
0.80x Mag.



83D biconcave lens facilitates viewing the fundus in an air filled vitreous cavity in phakic and pseudophakic eyes.




VMD

Magnifying lens
30° static F.O.V.
1.49x Mag.



For detailed examination and minute surgical manipulation of retinal membranes and pseudophakic eyes.




V30PD

30° Prism lens
33° static F.O.V.
1.02x Mag.



Provides visualization of the posterior peripheral fundus and vitreous beyond the equator with minimal distortion in phakic, aphakic and pseudophakic eyes.



Single Use

Volk Lenses

YAG Laser Lenses

VCAPSD10

Capsulotomy Laser Lens
Box of 10



VCAPSD10



VIRIDD10

Images not to scale

VIRIDD10

Iridectomy Laser Lens
Box of 10

Gonio Lenses

V3MIRD10

Three Mirror Gonio Lens
Box of 10



V3MIRD10



V4MIRD10

Images not to scale

V4MIRD10

Four Mirror Gonio Lens
Box of 10

SLT Lenses

VSLTD10

SLT Lens
Box of 10



VSLTD10

Images not to scale

Accessories

A7846

Precision Optical Lens Cleaner



Box of 24


Punctal Occlusion

Ultraplug Punctal Plugs

Silicone 'grommet' style punctum plugs pre-mounted on a combined applicator and dilator.

A10211



Punctal plugs, small, 0.4mm (silicone rubber) and applicator



Pack of 2


A10212

Punctal plugs, medium, 0.6mm (silicone rubber) and applicator



A10213

Punctal plugs, large, 0.8mm (silicone rubber) and applicator



A10225

Punctal gauging system, double ended 0.4mm/0.5mm, 0.6mm/0.7mm, 0.8mm/0.9mm, titanium



Set of 3
Reusable

Single Use

Punctal Occlusion

OPAQUE Herrick Lacrimal Plugs

For long term occlusion therapy. Mounted on disposable inserters. Radio and light opaque.

A10300

0.3mm diameter

A10302

0.5mm diameter

A10304

0.7mm diameter



Box of 2

Dissolvable OPAQUE Herrick Lacrimal Plugs

For medium term occlusion therapy. Dissolves in approx. 4-6 months. Radio and light opaque.

A10308

0.4mm diameter

A10306

0.5mm diameter



Box of 2

Accessories

A2080

Nettleship dilator



Reusable

A5908

Jewellers forceps,
No.5 110mm



Reusable

M9276

26G x 1.25" Curved
lacrimal cannula
(0.45mm x 32mm)



Box of 10

Blades and Knives

A1702

Scalpel handle, No.3



Reusable

A1704

No. 10 scalpel blades



Box of 100

A1706

No. 11 scalpel blades



A1710

No. 15 scalpel blades



A1725


Metal handle,
for Single Use blades



Reusable

A10130


Blade and tip,
15°



Box of 12


A10134

Blade and tip,
30°



A10136


Blade and tip,
45°



Blades and Knives

A10150


Sclerotome blade,
Type 57



Box of 6


A10154

Sclerotome blade,
angled Type 66



A10169


Myringotomy blade



Images not to scale

A10080


20G V-Lance



Box of 6


A10157

Sclerotome knife,
Type 57





A10182

Crescent knife,
angled,
bevel-up



A10169DK

Myringotomy knife



Images not to scale

Micro Keratome Knives, Angled

Box of 6


MSL12

1.20mm blade,
angled



MSL14

1.40mm blade,
angled



MSL16

1.60mm blade,
angled



MSL18

1.80mm blade,
angled



MSL20

2.00mm blade,
angled



MSL22

2.20mm blade,
angled




MSL23

2.30mm blade,
angled



MSL24

2.40mm blade,
angled



MSL25

2.50mm blade,
angled



Safety Knife version available on many models please suffix code with SK eg MSL25SK


Images not to scale

Micro Keratome Knives, Angled

Box of 6

MSL26

2.65mm blade,
angled



MSL27

2.75mm blade,
angled




MSL28

2.80mm blade,
angled



MSL29

2.90mm blade,
angled



MSL30

3.00mm blade,
angled



MSL32

3.20mm blade,
angled



MSL35


3.50mm blade,
angled



Short-Head Slit Knife

MSL**SH

Angled blade with
laser-etched 1.0mm
and 2.0mm width
incisions marks



Available in; 2.0mm, 2.1mm, 2.2mm, 2.3mm, 2.4mm, 2.5mm, 2.65mm, 2.75mm, 2.8mm, 2.9mm and 3.0mm (Catalogue number example MSL22SH)


Images not to scale

Safety Knife version available on many models please suffix code with SK eg MSL25SK

Crescent Knives

MCU26


2.30mm, bevel up,
angled



Box of 6


MCD26

2.30mm, bevel
down, angled



MCS26


2.30mm, straight



V-Lances


MVR19

19 gauge




MVR20

20 gauge




MVR20A

20 gauge,
angled




MVR23

23 gauge




MVR23A

23 gauge,
angled



MVR25

25 gauge



Images not to scale

Box of 6

LRI knives (Set depth blade)


MDP30 300 micron depth blade

MDP35 350 micron depth blade

MLR50 500 micron depth blade

MLR55 550 micron depth blade

MLR60 600 micron depth blade



Images not to scale


Blades and Knives



Round Slit Knives, Angled

MRS32


3.20mm blade,
angled



Box of 6


MRS35

3.50mm blade,
angled




MRS37

3.75mm blade,
angled





MRS41

4.10mm blade,
angled



MRS52

5.20mm blade,
angled



Microknives

Images not to scale

MST15


15° blade



Box of 6


MST15MK

15° blade with
laser-etched 1.00mm
and 1.20mm width
incisions marks



MST22

22.5° blade




MST30

30° blade



MST45

45° blade



Safety Knife version available on many models please suffix code with SK eg MST15SK

Images not to scale



Cellulose and PVA Products



Images not to scale

40-400

PVA spear
Box of 90
(18 x 5)

40-405


PVA spear
Box of 180
(36 x 5)

40-410

Cellulose spear
Box of 100
(20 x 5)

40-415

Cellulose spear
Box of 200
(20 x 10)



40-401

PVA point
Box of 90
(18 x 5)

40-406

PVA point
Box of 180
(18 x 10)

40-411

Cellulose point
Box of 100
(20 x 5)

40-416

Cellulose point
Box of 200
(20 x 10)

Images not to scale




40-420

Corneal light shield,
7mm diameter
Box of 20



40-431

Eye wick
Box of 20



40-430

Eye drain,
80cc capacity
Box of 20

40-435


Eye drain,
400cc capacity
Box of 10

Cellulose and PVA Products

40-462

Diamond knife
cleaning block
Box of 10

Images not to scale



40-900

Instrument wipe,
75mm x 75mm
Box of 20

40-820

LASIK shield,
9mm x 7mm
Box of 20

Images not to scale



40-822

LASIK shield,
8mm x 4mm
Box of 20

40-890

LASIK wipe
Box of 20


40-830

LASIK drain
Box of 20

A10002


Cellulose strip,
7mm x 35mm
Box of 500
(10 x 50)

Images not to scale




A10014

Cellulose point
Box of 500
(10 x 50)




A10016

Cellulose strip,
7mm x 35mm
Box of 150
(15 x 10)



A10000

Altomed Cellulose
spear
Box of 200
(20 x 10)



A10010

Altomed Cellulose
points
Box of 250
(25 x 10)





Miscellaneous

A2271

Iris retractor set

Set of 5



A2273

Iris retractor set

Set of 4

A10060

Eye shields,
clear, universal,
blister pack


Box of 100



C4100

Cautery,
low temp
(1300°f)

Box of 10



C4200

Cautery,
high temp
(2200°f)

Box of 10



Images not to scale

Miscellaneous

Ahmed Glaucoma Valves

The Ahmed Glaucoma Valve consists of a very small plate with a unique valve system which regulates eye pressure. The plate is attached to a tube which is inserted into the eye to drain excess fluid, thus reducing eye pressure. The valve adjusts the flow according to the pressure within the eye, ensuring there is no excessive drainage. The Ahmed Glaucoma Valve has been successfully demonstrated to work in all types of glaucoma cases, particularly neovascular, congenital and uveitic.

AGV-FP7

Flexible plate

AGV-FP8

Flexible plate small size

AGV-FX1

Flexible Bi-plate

AGV-S2


Standard size

TE


Tube extender

TI
Tube insertion forceps
(reusable)
AGV-FP7

(silicone)


AGV-S2

(polypropylene)



Images not to scale

Capsule Stain

A10330
Oftalblu
1ml vials
Box of 25Oftalblu intraocular staining
solution 0.1% Trypan blue.It is supplied sterile with a
perforation cap and security
seal.

Eyelid Implants / Weights

SGP06-SG28
Gold eyelid implants,
0.6 - 2.8g
in 0.2g increments

Sterile

Images not to scale

NST100
Tantalum sizing set
(0.6 - 1.8g)
Set of 7





Reusable

Labyrinthian lid loads are
thin, precisely tapered,
contoured 24 carat gold
weights that assist patients
with Lagophthalmos.The lid load adds sufficient
weight to allow gravity to
close the upper lid when
implanted.

Refractive Cannulas / Scleral Markers

Refractive Cannulas

Products designed for irrigating the cornea and stroma during laser refractive procedures.

Code	Gauge	Size	Description	Tip
				
M2804	25G x 7/8"	0.5 x 22mm	Irrigating Lasik Cannula. Formed 7mm from tip. Six 0.15mm ports at 2, 4 and 6mm from tip	
				
M2817	27G x 7/8"	0.4 x 22mm	Irrigating Lasik Cannula. Formed 9mm from tip. Four 0.15mm ports at 2 and 4mm from tip	
				
M2822	25G x 7/8"	0.5 x 22mm	Irrigating Lasik Cannula. Formed 7mm from tip. Six 0.15mm ports at 2, 4 and 6mm from closed tip	
				
M2824	26G x 7/8"	0.45 x 22mm	Irrigating Lasik Cannula. Formed 9.5mm from spatulated tip (Slade)	
				
M2801	25G x 7/8"	0.5 x 22mm	Irrigating Lasik Cannula. Flattened tip. Angled with two 0.25mm ports (Lin)	
				
M2802	25G x 1 1/8"	0.5 x 29mm	Irrigating Lasik Cannula. Formed with two 0.2mm side ports (Buratto)	
				
M3914	23G x 7/8"	0.64 x 22mm	Irrigating Lasik Cannula. Formed 10.5mm from flattened tip.	


Scleral Markers (Sterile and Non Sterile)

Code	Gauge	Size	Description	Tip
				
M7500	3.5mm/4.0mm	Scleral Marker. 80mm long	The scleral marker is designed to determine the injection or incision point through the pars plana, when located Scleral Marker. 80mm long against the limbus.	
				
M7510	3.5mm	Scleral Fixation Plate. 75mm long	The scleral fixation plate allows fixation of the globe and accurate injection through the pars plana when located against the limbus, using the 3.5mm slot.	

Anaesthetic Needles and Cannulas

Anaesthetic Cannulas

For injecting anaesthetic between Tenon's capsule and the sclera.

Code	Gauge	Size	Description	Tip
				
M0278	19G x 1"	1.1 x 25mm	Sub Tenon Cannula. Curved with flattened tip	
M0279	20G x 1"	0.9 x 25mm	Sub Tenon Cannula. Curved with flattened tip	

Anaesthetic Needles

For injecting anaesthetic around the globe.

Code	Gauge	Size	Description	Tip
				
M0275	 25G x 1 1/2"	0.5 x 38mm	Retrobulbar. Tri-facet	
M0637	 25G x 1 1/2"	0.5 x 38mm	Retrobulbar (Atkinson)	
M0638	 23G x 1 1/2"	0.64 x 38mm	Retrobulbar (Atkinson)	
				
M0641A	27G x 1 1/4"	0.4 x 32mm	Peribulbar (Atkinson)	
M0642	25G x 1 1/4"	0.5 x 32mm	Peribulbar (Atkinson)	
M0643	23G x 1 1/4"	0.64 x 32mm	Peribulbar (Atkinson)	
M0641	27G x 7/8"	0.4 x 22mm	Peribulbar (Atkinson)	
M0642A	25G x 7/8"	0.5 x 22mm	Peribulbar (Atkinson)	



Anterior Chamber Maintenance Cannulas

Anterior Chamber Maintainers

Self-retaining cannulas designed for the continual infusion of fluid to maintain the AC during cataract surgery or anterior vitrectomy.

Code	Gauge	Size	Description	Tip
				
M2658	23G x 9/64"	0.64 x 3.5mm	Anterior Chamber Maintainer. 150mm silicone tubing. Annular Rings allow cannula to be self retaining (Lewicky)	
M2657	20G x 9/64"	0.9 x 3.5mm	Anterior Chamber Maintainer. 150mm silicone tubing. Annular Rings allow cannula to be self retaining (Lewicky)	


Anterior Chamber Maintenance Cannulas

For maintaining and forming the anterior chamber during surgery by injecting or removing viscoelastic materials, BSS or air.

Code	Gauge	Size	Description	Tip
				
M2639	25G x 1"	0.5 x 25mm	Anterior Chamber Cannula. Angled 35°, 7.5mm from tip (Rycroft)	
M2270	30G x 7/8"	0.3 x 22mm	Anterior Chamber Cannula. Angled 35°, 4mm from tip (Rycroft)	
M2270A	30G x 7/8"	0.3 x 22mm	Anterior Chamber Cannula. Angled 35°, 8mm from tip (Rycroft)	
M2273	30G x 7/8"	0.3 x 22mm	Anterior Chamber Cannula. Angled 35°, 5.5mm from tip (Rycroft)	
M2273A	27G x 7/8"	0.4 x 22mm	Anterior Chamber Cannula. Angled 35°, 4.5mm from tip (Rycroft)	
M2279	27G x 7/8"	0.4 x 22mm	Anterior Chamber Cannula. Angled 45°, 11mm from tip (Rycroft)	
M2275	25G x 7/8"	0.5 x 22mm	Anterior Chamber Cannula. Angled 35°, 5mm from tip (Rycroft)	
M2273G	23G x 7/8"	0.64 x 22mm	Anterior Chamber Cannula. Angled 35°, 5.5mm from tip (Rycroft)	
M2642	21G x 7/8"	0.81 x 22mm	Anterior Chamber Cannula. Angled 35°, 7.5mm from tip (Rycroft)	
M2273B	20G x 7/8"	0.9 x 22mm	Anterior Chamber Cannula. Angled 45°, 7.5mm from tip (Rycroft)	


Transzonular Injection Cannula

For the intraoperative injection of steroid-antibiotic combinations.

Code	Gauge	Size	Description	Tip
				
M2280	27G x 7/8"	0.4 x 22mm	Transzonular Injection Cannula. Angled 33°, 2mm from tip.	




Wash Out Cannulas / Viscoelastic Cannulas

Wash Out Cannulas

Code	Gauge	Size	Description	Tip
				
M9280	19G x 1"	1.1 x 25mm	Washout Cannula. Angled 32°, 8mm from tip (Bishop Harmon). Tip flattened in the horizontal plane	
M9281	16G x 1"	1.6 x 25mm	Washout Cannula. Angled 30°, 6.5mm from the tip. Tip flattened in the horizontal plane	




Viscoelastic Cannulas


For injection of viscoelastic products into the anterior or posterior chamber.

Code	Gauge	Size	Description	Tip
				
M3926A	 27G x 7/8"	0.4 x 22mm	Viscoelastic Cannula. Angled at 35°, 9mm from tip	
M3927	 27G x 7/8"	0.4 x 22mm	Viscoelastic Cannula. Angled at 25°, 8mm from tip	
M2272	 26G x 7/8"	0.45 x 22mm	Viscoelastic Cannula. Angled at 35°, 9mm from tip	
M2639A	 25G x 7/8"	0.5 x 22mm	Viscoelastic Cannula. Angled at 35°, 9mm from tip	
M3909	 23G x 7/8"	0.64 x 22mm	Viscoelastic Cannula. Angled at 45°, 9mm from tip	
M3910A	 21G x 7/8"	0.81 x 22mm	Viscoelastic Cannula. Angled at 45°, 9mm from tip	
				
M3910	 21G x 7/8"	0.81 x 22mm	Viscoelastic Cannula. Straight	
				
M3917	 27G x 7/8"	0.4 x 22mm	Viscoelastic Cannula. Angled and formed 11mm from tip	



Irrigating cystotomes are designed for capsularhexis, intercapsular/endocapsular and can opener capsulotomy techniques. Formed cystotomes are shaped to conform to the lens. Short curved cystotomes are more suited to deep set eyes or small pupils. Capsularhexis cystotomes have a special tip which allows greater control of the flap when performing continuous curvilinear capsulotomies.

Formed Irrigating Cystotomes


Code	Gauge	Size	Description	Tip
				
M3630	 30G x 58"	0.3 x 16mm	Cystotome. Formed	
M3610B	 27G x 58"	0.4 x 16mm	Cystotome. Formed	
M3610	 25G x 58"	0.5 x 16mm	Cystotome. Formed	
M3610S	 25G x 58"	0.5 x 16mm	Cystotome. Short tip. Formed	
M3602	 25G x 58"	0.5 x 16mm	Cystotome. Side cutting. Formed	
				
M3610C	 25G x 58"	0.5 x 16mm	Cystotome. Tip angled 35°. Formed	
				
M3627A	 27G x 58"	0.4 x 16mm	Cystotome. Short curve. Formed.	
M3627	 25G x 58"	0.5 x 16mm	Cystotome. Short curve. Formed.	
				
M3626	 27G x 58"	0.4 x 16mm	Cystotome. Angled	
M3626A	 25G x 58"	0.5 x 16mm	Cystotome. Angled	
				
M3629	 25G x 58"	0.5 x 16mm	Cystotome. Reverse tip. Formed	
				
M3615	 27G x 78"	0.4 x 22mm	Long Cystotome. Formed.	



Formed Capsularhexis Cystotomes



Code	Gauge	Size	Description	Tip
				
M3000F	 27G x 58"	0.4 x 16mm	Cystotome. Formed	
M3000G	 25G x 58"	0.5 x 16mm	Cystotome. Formed	
				
M3000E	 27G x 58"	0.4 x 16mm	Cystotome. Short curve. Formed	
M3000D	 25G x 58"	0.5 x 16mm	Cystotome. Short curve. Formed	

Hydrodissection Cannulas






Hydrodissection Cannulas

A range of cannulas with tips designed to deliver fluids to aid the separation of the cortex from the nucleus and capsule.


Code	Gauge	Size	Description	Tip
				
M3919	27G x 7/8"	0.4 x 22mm	Hydrodissector. Angled 2mm tapered tip (Akahoshi)	
				
M2273C	27G x 7/8"	0.4 x 22mm	Hydrodissector. Curved. Tip flattened in vertical plane (Sauter)	
M2273D	25G x 7/8"	0.5 x 22mm	Hydrodissector. Curved. Tip flattened in vertical plane (Sauter)	
				
M2273E	27G x 7/8"	0.4 x 22mm	Hydrodissector. Curved. Tip flattened in horizontal plane (Helsinki)	
M2273F	25G x 7/8"	0.5 x 22mm	Hydrodissector. Curved. Tip flattened in horizontal plane (Helsinki)	
				
M2274	25G x 7/8"	0.5 x 22mm	Hydrodissector. Curved. Tip flattened in vertical plane (Kellan)	
				
M3992	30G x 7/8"	0.3 x 22mm	Hydrodissector. Angled 7mm from tip. Tip tapered and flattened in the horizontal plane	
M3902	27G x 7/8"	0.4 x 22mm	Hydrodissector. Angled 9mm from tip. Tip tapered and flattened in the horizontal plane	
M3900	25G x 7/8"	0.5 x 22mm	Hydrodissector. Angled 7mm from tip. Tip tapered and flattened in the horizontal plane	
				
M2273H	27G x 7/8"	0.4 x 22mm	Hydrodissector. Curved. Rounded tip	

Hydrodissection / Hydrodelineation / Viscoexpression Cannulas


Hydrodissector Cannulas

Code	Gauge	Size	Description	Tip
				
M3904	27G x 3/4"	0.4 x 19mm	Hydrodissector. Tip angled 90 degrees and flattened in the vertical plane (Chang)	
				
M3912	23G x 7/8"	0.64 x 22mm	Hydrodissector. Angled 30° 8mm from tip. Flattened 4mm from tip. Tip flattened in horizontal plane (Mackool)	
				
M3906A	25G x 7/8"	0.5 x 22mm	Hydrodissector. Round tip, 1.75mm wide. J-shaped tip	


Cortical Cleaving Cannulas

Code	Gauge	Size	Description	Tip
				
M3908	25G x 7/8"	0.5 x 22mm	Cortical Cleaving Hydrodissector. Tapered to 0.3mm at tip. Angled at 45°, 1mm from tip	
				
M3911	27G x 7/8"	0.4 x 22mm	Cortical Cleaving Hydrodissector. Angled at 45°, 11mm from tip. Flattened in horizontal plane 2.75mm from tip	
M3911A	27G x 7/8"	0.4 x 22mm	Cortical Cleaving Hydrodissector. Angled 45°, 11mm from round tip	

Hydrodelineation Cannulas

Code	Gauge	Size	Description	Tip
				
M3907	25G x 7/8"	0.5 x 22mm	Hydrodelineator. Tapered to 0.3mm at tip. Angled at 45°, 4.5mm from tip	

Viscoexpression Cannulas

Code	Gauge	Size	Description	Tip
				
M3899	25G x 7/8"	0.5 x 22mm	Viscoexpression Cannula. Flat tip angled 2mm from tip (Corydon)	
M3899A	25G x 7/8"	0.5 x 22mm	Viscoexpression Cannula. Round tip angled 1.8mm from tip (Corydon)	

Nucleus / Cortex Removal Cannulas


Nucleus Removal Cannulas

Irrigating Vectis are designed for irrigating and removing the nucleus and the larger pieces of cortical debris whilst maintaining the anterior chamber depth.

Code	Gauge	Size	Description	Tip
				
M4619	25G x 1 1/2"	0.5 x 38mm	Irrigating Vectis. 5mm wide, single port	
M4620	25G x 1 1/2"	0.5 x 38mm	Irrigating Vectis. 5mm wide, twin ports	

Cortex Removal Cannulas



For aspirating cortex debris when chamber is maintained with an anterior chamber maintainer or viscoelastic substance. The 'J' shaped cannulas are designed to retract the iris and aspirate difficult to reach cortex between the 10 to 2 O'clock position.

Code	Gauge	Size	Description	Tip
				
M6576	23G x 7/8"	0.64 x 22mm	Cortex Extractor. 0.3mm top port	
M6575	23G x 5/8"	0.64 x 16mm	Cortex Extractor. 0.3mm top port	
M6575A	23G x 5/8"	0.64 x 16mm	Cortex Extractor. 0.4mm top port	
M6576A	23G x 1/2"	0.64 x 12.5mm	Cortex Extractor. 0.3mm top port	
M6594	23G x 5/8"	0.64 x 16mm	Cortex Extractor. 0.3mm bevelled end port (Charleux)	
M6596	24G x 1/2"	0.56 x 12.5mm	Cortex Extractor. 0.3mm bevelled end port (Charleux)	
M6595	23G x 1/2"	0.64 x 12.5mm	Cortex Extractor. 0.3mm bevelled end port (Charleux)	



Capsule Polishing Cannulas

Capsule Polishing Cannulas



Designed for polishing and cleaning the residual cortex from the capsular bag. Silicone tipped cannulas provide gentle polishing and a traumatic insertion and removal. Micro-etched cannulas are slightly roughened for gentle polishing. Olive tipped cannulas allow polishing over a greater surface area with a reduced risk of posterior capsule tear.

Code	Gauge	Size	Description	Tip
				
M8604A	 23G x 7/8"	0.64 x 22mm	Capsule Polisher. Angled 9mm from tip. Aspiration through 0.3mm side port. Micro-etched tip	
M8604	 21G x 1"	0.81 x 25mm	Capsule Polisher. Angled 12.5mm from tip. Aspiration through 0.3mm side port. Micro-etched tip	
M8604B	 21G x 7/8"	0.81 x 22mm	Capsule Poilisher. Angled 9mm from tip. Aspiration through 0.3mm side port. Micro-etched tip.	
				
M8645	 23G x 7/8"	0.64 x 22mm	Capsule Scraper. Angled 0.5 & 7mm from tip. Micro-etched distal end.	

Silicone Tipped Polishers

Code	Gauge	Size	Description	Tip
				
M8605	 27G x 7/8"	0.4 x 22mm	Silicone tipped polisher. Angled	
M8605A	 23G x 7/8"	0.64 x 22mm	Silicone tipped polisher. Angled	
				
M8608	 27G x 7/8"	0.4 x 22mm	Silicone tipped polisher. Curved through 45°	

Olive Tipped Cannulas

Code	Gauge	Size	Description	Tip
				
M8622	25G x 1"	0.5 x 25mm	Olive tipped polisher. Angled 10mm from tip. Micro-etched tip.	
				
M8620C	23G x 1"	0.64 x 25mm	Olive tipped polisher. Curved	
M8621	23G x 1"	0.64 x 25mm	Olive tipped polisher. Curved. Micro-etched tip	

Irrigating/Aspirating Cannulas and Handpieces

Coaxial Cannulas and Handpieces

For the aspiration of cortex while simultaneously maintaining the AC with a continuous irrigating flow.

Code	Gauge	Size	Description	Tip
				
M6634	18G/21G	1.24/0.81mm	Coaxial I & A Cannula. 0.3mm side aspiration port on inner cannula	
				
M6611	21G x 1 1/2"	0.81 x 38mm	I/A Straight Cannula	
				
M5590	17/24G	1.47/0.56mm	I/A Handpiece for 2.2mm incisions. Straight tip	
M5592	17/24G	1.47/0.56mm	I/A Handpiece for 2.2mm incisions. 45° angled tip	
M5593	17/24G	1.47/0.56mm	I/A Handpiece for 2.2mm incisions. 60° angled tip	
M5594	17/24G	1.47/0.56mm	I/A Handpiece for 2.2mm incisions. 90° angled tip	
				
M5591	17/24G	1.47/0.56mm	I/A Handpiece for 2.2mm incisions. Curved tip	
				
M5580	16/21G	1.6/0.81mm	I/A Handpiece for 2.7/2.8mm incisions. Straight tip	
M5582	16/21G	1.6/0.81mm	I/A Handpiece for 2.7/2.8mm incisions. 45° angled tip	
M5583	16/21G	1.6/0.81mm	I/A Handpiece for 2.7/2.8mm incisions. 60° angled tip	
M5584	16/21G	1.6/0.81mm	I/A Handpiece for 2.7/2.8mm incisions. 90° angled tip	
				
M5581	16/21G	1.6/0.81mm	I/A Handpiece for 2.7/2.8mm incisions. Curved tip.	
				
M5560	21G	Soft Outer	I/A Handpiece for 2.7/2.8mm incisions. Straight tip	
M5561	21G	Soft Outer	I/A Handpiece for 2.7/2.8mm incisions. 45° angled tip	
M5562	21G	Soft Outer	I/A Handpiece for 2.7/2.8mm incisions. 60° angled tip	
M5563	21G	Soft Outer	I/A Handpiece for 2.7/2.8mm incisions. 90° angled tip	

Irrigating/Aspirating Cannulas and Handpieces

Code	Gauge	Size	Description	Tip
				
M5570	24G	Soft Outer	I/A Handpiece for 2.2mm incisions. Straight tip	
M5571	24G	Soft Outer	I/A Handpiece for 2.2mm incisions. 45° angled tip	
M5572	24G	Soft Outer	I/A Handpiece for 2.2mm incisions. 60° angled tip	
M5573	24G	Soft Outer	I/A Handpiece for 2.2mm incisions. 90° angled tip	

Simcoe Irrigating / Aspirating Cannulas

For removing cortical debris or viscoelastic substances while maintaining the anterior chamber during extracapsular and intercapsular/endocapsular cataract procedures.

Code	Gauge	Size	Description	Tip
				
M6600	● 23G/23G	0.64/0.64mm	I/A Cannula Standard. 0.4mm port. 200mm silicone tubing 2.65 OD x 0.65 ID (Simcoe)	
M6602	● 23G/23G	0.64/0.64mm	I/A Cannula Standard. 0.4mm port. 270mm silicone tubing 2.65 OD x 0.65 ID (Simcoe)	
				
M6605	● 23G/23G	0.64/0.64mm	I/A Cannula Standard. 0.3mm port. 200mm silicone tubing 2 OD x 1 ID (Simcoe)	

Simcoe Irrigating / Aspirating Cannulas

Code	Description
	
M1000	Irrigating Microsurgical Handle with female luer taper and a male luer lock connector

Irrigating/Aspirating Cannulas and Handpieces

Bimanual I/A Handpieces

All Bimanual handles have a coloured indicator to show the gauge size in accordance with ISO 6009.

Code	Gauge	Size	Description	Tip
------	-------	------	-------------	-----



M5501	23G x 5/8"	0.64 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.4mm side port and micro-etched tip	
M5502	23G x 5/8"	0.64 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port and micro-etched tip	
M5512	21G x 5/8"	0.81 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port and micro-etched tip	
M5520	21G x 5/8"	0.81 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.4mm side port and micro-etched tip	
M5511	20G x 3/4"	0.9 x 19mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port and micro-etched tip	
M5518	23G x 5/8"	0.64 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port	
M5515	21G x 7/8"	0.81 x 22mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port	



M5519	23G x 5/8"	0.64 x 16mm	Bimanual Irrigation Handpiece. Ultra Thin Wall. Curved with 45° bevelled tip	
M5505	21G x 5/8"	0.81 x 16mm	Bimanual Irrigation Handpiece. Curved with 45° bevelled tip	
M5513	23G x 5/8"	0.64 x 16mm	Bimanual Irrigation Handpiece. Thin Wall. Curved with twin 0.4mm side ports	
M5521	21G x 5/8"	0.81 x 16mm	Bimanual Irrigation Handpiece. Thin Wall. Curved with twin 0.5mm side ports	
M5514	21G x 5/8"	0.81 x 16mm	Bimanual Irrigation Handpiece. Curved with twin 0.5mm side ports	
M5510	20G x 3/4"	0.9 x 19mm	Bimanual Irrigation Handpiece. Curved with twin 0.5mm side ports	
M5523	23G x 5/8"	0.64 x 16mm	Bimanual Irrigation Handpiece. Ultra Thin Wall. Curved with open end	

Bimanual I/A Systems



A combination set of one Irrigating Handle and one Aspirating Handle manufactured using thin wall stainless steel tubing.

Code	Gauge	Size	Description	Tip
------	-------	------	-------------	-----





M5537 - set	20G x 3/4"	0.9 x 19mm	Bimanual Irrigation Handpiece. Curved with twin 0.5mm side ports	
	21G x 7/8"	0.81 x 22mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port	

Bimanual I/A Systems

Code	Gauge	Size	Description	Tip
				
M5531 - set	21G x 58"	0.81 x 16mm	Bimanual Irrigation Handpiece. Thin Wall. Curved with twin 0.5mm side ports	
	21G x 58"	0.81 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.4mm side port and micro-etched tip	
				
M5331 - set	21G x 58"	0.81 x 16mm	Bimanual Irrigation Handpiece. Curved with 2 x 0.4mm irrigation ports and tapered tip	
	21G x 58"	0.81 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm aspiration port and micro-etched tapered tip	
				
M5535 - set	23G x 58"	0.64 x 16mm	Bimanual Irrigation Handpiece. Thin Wall. Curved with twin 0.4mm side ports	
	23G x 58"	0.64 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port and micro-etched tip	
M5536 - set	23G x 58"	0.64 x 16mm	Bimanual Irrigation Handpiece. Ultra Thin Wall. Curved with 45° bevelled tip	
	23G x 58"	0.64 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port	
M5538 - set	25G x 58"	0.5 x 16mm	Bimanual Irrigation Handpiece. Ultra Thin Wall. Curved with 45° bevelled tip	
	25G x 58"	0.5 x 16mm	Bimanual Aspiration Handpiece. Curved with 0.3mm side port	



Bimanual Cannulas

Code	Gauge	Size	Description	Tip
				
M6524	21G x 58"	0.81 x 16mm	Bimanual Aspiration Cannula. Curved with 0.3mm top port	
				
M5405	21G x 58"	0.81 x 16mm	Bimanual Irrigation Cannula. Curved with 45° bevelled tip	

Lacrimal Cannulas and Dacryocystorhinostomy (DCR) Sets

Lacrimal Cannulas


Cannulas for probing and irrigating the lacrimal ducts.

Code	Gauge	Size	Description	Tip
				
M9275	 26G x 3/4"	0.45 x 20mm	Lacrimal Cannula. Short straight	
				
M9276	 26G x 1 1/4"	0.45 x 32mm	Lacrimal Cannula. Curved	
M9277	 25G x 1 1/8"	0.5 x 29mm	Lacrimal Cannula. Curved	
M9274	 21G x 1 1/2"	0.81 x 38mm	Lacrimal Cannula. Curved	
				
M9278	 25G x 7/8"	0.5 x 22mm	Lacrimal Cannula. Straight with 0.2mm side port	
				
M9272	 23G x 1/2"	0.64 x 13mm	Lacrimal Cannula. Straight with 0.3mm side port. Tapered tip.	

Dacryocystorhinostomy (DCR) Sets

Sterimedix DCR sets offer the latest manufacturing technology and design. Improved attachment of the silicone rod to the malleable probes offers the practitioner the best overall device for DCR surgery. Tapered tips facilitate entry into the lacrimal ducts.

Code	Gauge	Size	Description	Tip
M9290	20G x 1 3/4"	0.9 x 45mm	Straight DCR Set with 390mm silicone rod	
M9291	20G x 1 3/4"	0.9 x 45mm	Angled DCR Set with 390mm silicone rod	



Vitreoretinal Extension Tube & Self Retaining Infusion Cannulas


High Pressure Extension Tube

Code	Size	Description
		
M4463	125mm	High Pressure Extension Tube

Self Retaining Infusion Cannulas

For simple insertion and removal without the need for a retaining suture. Designed for use directly into the eye without the need for a trocar.

Code	Gauge	Size	Description	Tip
				
M4464	20G	4.3mm 300mm overall length	Self Retaining Infusion Cannula. Angled	
				
M4467A	20G	6mm 150mm overall length	Self Retaining Infusion Cannula. Straight	




One-step Trocar Sets

Code	Gauge	Description
------	-------	-------------



M4510	23G	One-step Trocar Set containing: 3 x self-sealing valved trocars 1 x 3.5mm scleral fixation plate 1 x 185mm BSS infusion line
M4515	25G	One-step Trocar Set containing: 3 x self-sealing valved trocars 1 x 3.5mm scleral fixation plate 1 x 185mm BSS infusion line




Vitreoretinal Cannulas


Vitreoretinal Cannulas

Code	Gauge	Size	Description	Tip
				
M4477	 25G x 1 1/4"	0.5 x 32mm	Aspiration Cannula	
M4478	 23G x 1 1/4"	0.64 x 32mm	Aspiration Cannula	
M4478B	 23G UTW x 1 1/4"	0.64 x 32mm	Aspiration Cannula. Ultra Thin Wall	
M4471	 20G x 1 1/4"	0.9 x 32mm	Aspiration Cannula	
				
M4410	 20G	0.9mm	Subretinal Fluid Cannula. Straight. Silicone sleeve extends 6mm from end of cannula	

Viscous Fluid Injection Cannulas

Code	Gauge	Size	Description	Tip
				
M4454	 20G x 1/4"	0.90 x 7mm	Viscous Fluid Injection Cannula	
M4453	 25G x 1/4"	0.5 x 7mm	Viscous Fluid Injection Cannula	
M4452	 23G x 1/4"	0.64 x 8mm	Viscous Fluid Injection Cannula. Bevelled tip	
				
M4447	 23G x 1/4"	0.64 x 4mm	Viscous Fluid Injection Cannula	
M4448	 25G x 1/4"	0.5 x 4mm	Viscous Fluid Injection Cannula	

Dual-Bore Cannulas

Code	Gauge	Size	Description	Tip
				
M4474	20G/25G	0.9/0.5mm	Perfluorocarbon Coaxial (Dual Bore) Cannula	
M4499	23G/30G	0.64/0.3mm	Perfluorocarbon Coaxial (Dual Bore) Cannula	
M4499A	25G/30G	0.5/0.3mm	Perfluorocarbon Coaxial (Dual Bore) Cannula	

Vitreoretinal Cannulas and Vitreoretinal Handles

Silicone Brush Cannulas

Code	Gauge	Size	Description	Tip
				
M4449	25G x 1 3/8"	0.5 x 35mm	Silicone Brush. Silicone extends 0.7mm beyond end of cannula	
M4459	23G x 1 1/4"	0.64 x 32mm	Silicone Brush. Silicone extends 0.7mm beyond end of cannula	
M4475	20G x 1 1/4"	0.9 x 32mm	Silicone Brush. Silicone extends 1.25mm beyond end of cannula	
M4475A	20G x 1 1/4"	0.9 x 32mm	Silicone Brush. Silicone extends 0.6mm beyond end of cannula	

Soft Tip Cannulas

Code	Gauge	Size	Description	Tip
				
M4456	25G x 1 1/4"	0.5 x 32mm	Soft Tip. Silicone extends 1mm beyond end of cannula	
M4458	23G x 1 1/4"	0.64 x 32mm	Soft Tip. Silicone extends 2mm beyond end of cannula	
M4457	20G x 1 1/4"	0.9 x 32mm	Soft Tip. Silicone extends 1mm beyond end of cannula	

Backflush Handles and Combination Sets

Code	Gauge	Size	Description	Tip
				
M4489A	M4489A 'Active/Passive' Backflush Handle with finger control. Allows conversion from active to passive aspiration.			
				
M4504 - set	 25G x 1 1/4"	0.5 x 32mm	'Active/Passive' Backflush Handle with 25G Soft Tip Cannula	
M4530 - set	 23G x 1 1/4"	0.64 x 32mm	'Active/Passive' Backflush Handle with 23G Soft Tip Cannula	
M4531 - set	 20G x 1 1/4"	0.9 x 32mm	'Active/Passive' Backflush Handle with 20G Soft Tip Cannula	

Vitreoretinal Heavy Liquid Infusion Handles

Code	Description
	
M4498	Heavy Liquid Infusion Handle with finger control and a luer lok connection

Backflush Handles and Combination Sets

Code	Gauge	Size	Description	Tip
				
M4489	Backflush Handle for Passive Aspiration with finger control and a luer lok cannula connection			
				
				Soft Tip
M4501 - set	 20G x 1 1/4"	0.9 x 32mm	Backflush Handle for Passive Aspiration with 20G Soft Tip Cannula set	
M4502 - set	 23G x 1 1/4"	0.64 x 32mm	Backflush Handle for Passive Aspiration with 23G Soft Tip Cannula set	
M4503 - set	 25G x 1 1/4"	0.5 x 32mm	Backflush Handle for Passive Aspiration with 25G Soft Tip Cannula set	
				
				Silicone Brush
M4491 - set	 20G x 1 1/4"	0.9 x 32mm	Backflush Handle for Passive Aspiration with 20G Silicone Brush set	
M4492 - set	 23G x 1 1/4"	0.64 x 32mm	Backflush Handle for Passive Aspiration with 23G Silicone Brush set	
				
				Aspiration Cannula
M4496 - set	 20G x 1 1/4"	0.9 x 32mm	Backflush Handle for Passive Aspiration with 20G Aspiration Cannula set	
M4497 - set	 23G x 1 1/4"	0.64 x 32mm	Backflush Handle for Passive Aspiration with 23G Aspiration Cannula set	
M4500 - set	 25G x 1 1/4"	0.5 x 32mm	Backflush Handle for Passive Aspiration with 25G Aspiration Cannula set	

Backflush Handles and Combination Sets


Code	Gauge	Size	Description	Tip
				
M4493			'Charles' Flute Handle	
				
M4494 - set	20G x 1 1/4"	0.9 x 32mm	'Charles' Flute Handle with 20G Aspiration Cannula set	
M4495 - set	23G x 1 1/4"	0.64 x 32mm	'Charles' Flute Handle with 23G Aspiration Cannula set	

Drapes and Fluid Collection

Sterimedix drapes are offered in various sizes and materials and are suitable for the left or right eye. They are either waterproof or moisture resistant and the incise offers sufficient adhesive without being too aggressive on removal.


AMD 1020

40cm x 40cm
(16" x 16")
Aperture Drape
Flexible polymer
Incise Area:
110 x 110mm
Aperture:
65mm diameter
Box of 10




AMD 1021

40cm x 40cm
(16" x 16")
Incise Drape
Flexible polymer
Incise Area:
110 x 120mm
Box of 10




AMD 1061

122cm x 127cm
(48" x 50")
Incise Drape with
Fluid Collection
Pouch, flexible
polymer
Incise Area:
110 x 120mm
Box of 5




AMD 1062

122cm x 127cm
(48" x 50")
Incise Drape with
Fluid
Collection Pouch,
non-woven
Incise Area:
110 x 120mm
Box of 5



AMD 1063

122cm x 127cm
(48" x 50")
Aperture Drape
with Fluid Collection
Pouch, non-woven
Incise Area:
110 x 120mm
Aperture:
44 x 75mm
Box of 5



Colour Key

● Fluid collection pouch

● Adhesive area

● Flexible polymer

● Soft non-woven fabric

Index

A

Ahmed Glaucoma Valves.....	122
Altomed Needleholder.....	45
Anaesthetic Cannulas	124
Anaesthetic Needles	124
Anterior Chamber Maintainers.....	125
Anterior Chamber Maintenance Cannulas.....	125

B

Backflush Handles and Combination Sets	141, 142
Barraquer Needleholders	41, 42
Barron Marking Punches	98
Bimanual Cannulas	136
Bimanual I/A Handpieces	29, 30, 135
Bimanual I/A Systems	135, 136
Bipolar	57, 102
Bipolar Cables.....	59
Machine Fittings	59
Bipolar Forceps.....	58
Bipolar Pencils	58
Birks Needleholder	45
Blades and Knives.....	113, 114
Bonn Scissors	5
Bottle Parts	95

C

Calipers.....	54, 55
Callahan (Cox) Tube	52
Cannulas	
Anaesthetic Cannulas	124
Anterior Chamber Maintenance Cannulas.....	125
Bimanual Cannulas	136
Capsule Polishing Cannulas.....	132
Coaxial Cannulas and Handpieces	133
Cortex Removal Cannulas	131
Cortical Cleaving Cannulas	130
Dual-Bore Cannulas	140
Hydrodelineation Cannulas	130
Hydrodissection Cannulas	129
Hydrodissector Cannulas.....	130
Lacrimal Cannulas	137
Nucleus Removal Cannulas.....	131
Olive Tipped Cannulas	132
Refractive Cannulas	123
Self Retaining Infusion Cannulas	138
Silicone Brush Cannulas	141
Silicone Tipped Polishers	132
Simcoe Irrigating / Aspirating Cannulas	134
Soft Tip Cannulas.....	141
Transzonular Injection Cannula.....	125
Viscoelastic Cannulas	126
Viscoexpression Cannulas	130
Viscous Fluid Injection Cannulas	140
Vitreoretinal Cannulas	140

Wash Out Cannulas.....	126
Capsule Polishing Cannulas.....	132
Capsule Scissors	9
Capsule Stain	122
Capsulotomy	
Formed Capsularhexis Cystotomes	128
Formed Irrigating Cystotomes.....	127
Castroviejo Corneal Scissors.....	9
Castroviejo Micro Needleholders	39
Castroviejo Micro Scissors.....	10
Castroviejo Needleholders.....	43, 44
Castroviejo Scissors	7
Catford Micro Needleholders.....	40
CDCR	53
CDCR Tubes	52
Callahan (Cox) Tube	52
Double Collar Tube	52
Gladstone Putterman Tube	52
Lester Jones Tubes, Straight and Angled.....	52
Cellulose Products	119, 120
Circling Bands.....	105
Citelli Punch.....	46
Classic Indirect Bio Lenses.....	90, 108
Coaxial Cannulas and Handpieces.....	133
Conjunctival Scissors.....	6
Corneal Burrs.....	56
Corneal Markers	85, 86
Corneal Punches.....	87, 98, 99
Barron Marking Punches.....	98
Corneal Trephines	87, 98, 99
Hessburg-Barron Vacuum Trephines	98
Cortex Removal Cannulas	131
Cortical Cleaving Cannulas	130
Crawford Lacrimal Intubation System	51, 100
Crescent Knives.....	117

D

DCR (Dacryocystorhinostomy)	46 - 51, 100
Citelli Punch.....	46
Crawford Lacrimal Intubation System	51, 100
Dacryocystorhinostomy Sets	137
DCR Trephines.....	48
Hajeks Punch.....	46
Malhotra Endonasal Punch	47
O'Donoghue DCR Bodkins.....	100
Sphenoidal Punch	48
Depressors	54, 55
De Wecker Scissors.....	7
Diamond Knives.....	12 - 14
Digital Indirect Bio Lenses.....	90
Digital Slit Lamp Lenses	90
Double Collar Tube	52
Drapes	143
Dressing Scissors	5
Dual-Bore Cannulas	140

E

Enucleation.....	88
Enucleation Scissors.....	4
Eyelid Implants / Weights.....	122
Eye Scissors.....	1

F

Fluid Collection	143
Flute Handles	107
Forceps.....	60 - 80
Formed Capsularhexis Cystotomes	128
Formed Irrigating Cystotomes.....	127

G

Gladstone Putterman Tube.....	52
Global Fixation	85, 86
Gonio Lenses.....	92, 110
Grooved Scleral Sponge	104
Grooved Strips.....	105

H

Hajeks Punch.....	46
Herrick Lacrimal Plugs	112
Herrick Lacrimal Plugs Dissolvable	112
Hessburg-Barron Vacuum Trephines.....	98
High Pressure Extension Tube	138
Hooks, Probes and Chops.....	15 - 27
Hydrodelineation Cannulas.....	130
Hydrodissection Cannulas	129
Hydrodissector Cannulas.....	130

I

Indirect Contact Laser Lenses.....	91
Infusion Cannula.....	107
Intertip I/A System	28
Iris Scissors.....	2
Irrigation / Aspiration.....	28 - 30
Bimanual I/A Handpieces	29, 30, 135
Bimanual I/A Systems.....	135, 136
Coaxial Cannulas and Handpieces	133
Intertip I/A System	28
Simcoe Irrigating / Aspirating Cannulas	134

K

Katzin Corneal Scissors.....	10
------------------------------	----

L

Lacrimal Cannulas.....	137
Lester Jones Tubes	52
Lid Surgery.....	81 - 84
Lim Needleholders.....	45

LRI knives.....	117
Lucite Conformers.....	101

M

Machine Fittings	59
Malhotra Endonasal Punch.....	47
Mani Knives	115 - 118
Crescent Knives.....	117
LRI knives.....	117
Micro Keratome Knives.....	115, 116
Microknives.....	118
Round Slit Knives.....	118
Short-Head Slit Knife	116
Micro Keratome Knives	115, 116
Microknives	118
Microwash Trays.....	93
Sterilization Trays.....	93

N

Needleholders	38 - 45
Altomed Needleholder.....	45
Barraquer Needleholders	41, 42
Birks Needleholder.....	45
Castroviejo Micro Needleholders	39
Castroviejo Needleholders.....	43, 44
Catford Micro Needleholders.....	39, 40
Lim Needleholders.....	45
Silcocks Needleholder	45
Titanium Micro Needleholders.....	38
Troutman Needleholders.....	41
Nucleus Removal Cannulas	131

O

Oculoplastic Oncology	101
Lucite Conformers.....	101
Silicone Conformers	101
Silicone Spheres.....	101
O'Donoghue DCR Bodkins.....	100
Olive Tipped Cannulas.....	132
One-step Trocar Sets	139
Operating Stools.....	95
Operating Table Pillows.....	93
Orbit.....	88
Oval Scleral Sponge.....	103

Index

P

Part Thickness Scleral Sponge	104
Pillows.....	93
Operating Table Pillows.....	93
Polymer Sterilization Trays.....	94
Pooley Conjunctival Scissors.....	4
Pressure Reducers	56
Punches.....	54, 55
Punctal Occlusion.....	111, 112
Accessories.....	112
Herrick Lacrimal Plugs	112
Herrick Lacrimal Plugs Dissolvable	112
Ultraplug Punctal Plugs.....	111
PVA Products.....	119, 120

Q

Q-Rinse.....	95
Bottle Parts.....	95

R

Refractive Cannulas	123
Retractors and Rings	36, 37
Round Scleral Sponge	103
Round Slit Knives.....	118

S

Scholda Sponge.....	103
Scissors.....	1 - 10
Bonn Scissors	5
Capsule Scissors	9
Castroviejo Corneal Scissors.....	9
Castroviejo Micro Scissors	10
Castroviejo Scissors.....	7
Conjunctival Scissors.....	6
De Wecker Scissors.....	7
Dressing Scissors	5
Enucleation Scissors.....	4
Eye Scissors	1
Iris Scissors	2
Katzin Corneal Scissors.....	10
Pooley Conjunctival Scissors	4
Stevens Tenotomy Scissors	3, 4
Stitch Scissors	5
Strabismus Scissors.....	2, 3
Vannas Scissors.....	8
Werbs Scissors	4
Westcott Tenotomy Scissors.....	7
Williamson Noble Scissors.....	6
Scleral Markers	123
Scleral Sponge Accessories.....	104
Self Retaining Infusion Cannulas.....	138
Short-Head Slit Knife	116
Silcocks Needleholder	45

Silicone Brush Cannulas.....	141
Silicone Conformers	101
Silicone Lace.....	105
Silicone Spheres.....	101
Silicone Tipped Polishers	132
Simcoe Irrigating / Aspirating Cannulas	134
Sleeves.....	105
Slit Lamp Lenses	90
SLT Lenses	110
Soft Tip Cannulas.....	141
Specula.....	31 - 35
Sphenoidal Punch.....	48
Sterilization Trays.....	93, 94
Polymer Sterilization Trays	94
Stevens Tenotomy Scissors	3, 4
Stitch Scissors	5
Strabismus Scissors.....	2, 3
Surgeons Stools	95
Operating Stools.....	95
Surgical Knives and Blades	11

T

Tantalum Clips	105
Titanium Micro Needleholders.....	38
Transzonular Injection Cannula	125
Troutman Needleholders	41
Tunnel Scleral Sponge	103
Tyre Implants	106

U

Ultraplug Punctal Plugs	111
-------------------------------	-----

V

Vannas Scissors	8
Viscoelastic Cannulas.....	126
Viscoexpression Cannulas.....	130
Viscous Fluid Injection Cannulas	140
Vitreotomy Lenses	109
Vitreoretinal.....	89, 103 - 109
Backflush Handles and Combination Sets.....	141, 142
Circling Bands	105
Flute Handles	107
Grooved Scleral Sponge.....	104
Grooved Strips	105
High Pressure Extension Tube	138
Infusion Cannula.....	107
Miscellaneous Components.....	106
Oblong Scleral Sponge.....	103
One-step Trocar Sets	139
Oval Scleral Sponge.....	103
Part Thickness Scleral Sponge	104
Round Scleral Sponge.....	103
Scholda Sponge	103
Scleral Sponge Accessories.....	104

Silicone Lace.....	105
Sleeves.....	105
Tantalum Clips	105
Tunnel Scleral Sponge	103
Tyre Implants.....	106
Vitreotomy Lenses	109
Vitreoretinal Heavy Liquid Infusion Handles..	141
Vitreoretinal Cannulas	140
Vitreoretinal Heavy Liquid Infusion Handles....	141
Volk Lenses	90 - 92, 110
Accessories.....	92, 110
Classic Indirect Bio Lenses.....	90, 108
Digital Indirect Bio Lenses	90
Digital Slit Lamp Lenses	90
Gonio Lenses	92, 110
Indirect Contact Laser Lenses	91
Slit Lamp Lenses	90
SLT Lenses	110
YAG Laser Lenses	91, 110

W

Wash Out Cannulas.....	126
Werbs Scissors	4
Westcott Tenotomy Scissors.....	7
Williamson Noble Scissors.....	6

Y

YAG Laser Lenses	91, 110
------------------------	---------

Numerical Index

2058-71.....	59	51-903-8.75.....	99	A1145.....	3	A1650.....	11
2058-72.....	59	51-903-9.00.....	99	A1150.....	3	A1660.....	11
2063-71.....	59	51-903-9.25.....	99	A1155.....	3	A1662.....	11
2063-72.....	59	51-903-9.50.....	99	A1160.....	3	A1696.....	11
28-0174.....	51, 100	A10000.....	120	A1165.....	4	A1702.....	11, 113
28-0175.....	51, 100	A10002.....	120	A1170.....	4	A1704.....	113
28-0184.....	51, 100	A10010.....	120	A1180.....	4	A1706.....	113
28-0185.....	51, 100	A10014.....	120	A1185.....	4	A1710.....	113
28-0186.....	51, 100	A10016.....	120	A1190.....	5	A1725.....	11, 113
28-0189.....	51	A10060.....	121	A11908S.....	52	A2012.....	15
28-0190.....	51	A10080.....	114	A11930.....	52	A2018.....	15
40-400.....	119	A10130.....	113	A11932.....	52	A2020.....	15
40-401.....	119	A10134.....	113	A11934.....	52	A2022.....	15
40-405.....	119	A10136.....	113	A11936.....	52	A2024.....	15
40-406.....	119	A10150.....	114	A11937.....	52	A2026.....	15
40-410.....	119	A10154.....	114	A11937A.....	52	A2038.....	16
40-411.....	119	A10157.....	114	A11937B.....	52	A2040.....	16
40-415.....	119	A10169.....	114	A11940.....	53	A2042.....	16
40-416.....	119	A10169DK.....	114	A1195.....	5	A2044.....	16
40-420.....	119	A10182.....	114	A1235A.....	5	A2046.....	16
40-430.....	119	A10211.....	111	A1240.....	5	A2048.....	16
40-431.....	119	A10212.....	111	A1244.....	5	A2050.....	17
40-435.....	119	A10213.....	111	A1250.....	6	A2080.....	17, 112
40-462.....	120	A10225.....	111	A1251.....	6	A2087.....	17
40-820.....	120	A10300.....	112	A1252.....	6	A2091.....	17
40-822.....	120	A10302.....	112	A1300.....	6	A2094.....	17
40-830.....	120	A10304.....	112	A1310.....	6	A2110.....	18
40-890.....	120	A10306.....	112	A1312.....	6	A2122.....	18
40-900.....	120	A10308.....	112	A1314.....	6	A2135.....	18
51-900-6.00.....	99	A10330.....	122	A1315.....	7	A2142.....	18
51-900-6.25.....	99	A1060.....	1	A1325.....	7	A2144.....	18
51-900-6.50.....	99	A1065.....	1	A1350.....	7	A2145.....	18
51-900-6.75.....	99	A1070.....	1	A1355.....	7	A2152.....	19
51-900-7.00.....	99	A1075.....	1	A1360.....	7	A2154.....	19
51-900-7.25.....	99	A1100.....	1	A1365.....	7	A2158.....	19
51-900-7.50.....	99	A11020.....	95	A1380.....	8	A2160.....	19
51-900-7.75.....	99	A11022.....	95	A1385.....	8	A2162.....	19
51-900-8.00.....	99	A11024.....	95	A1395.....	8	A2170.....	19
51-900-8.25.....	99	A11026.....	95	A1400.....	8	A2184.....	19
51-900-8.50.....	99	A11028.....	95	A1405.....	8	A2196.....	20
51-900-8.75.....	99	A11029.....	95	A1415.....	8	A2198.....	20
51-900-9.00.....	99	A11044A.....	95	A1435.....	9	A2200.....	20
51-900-9.25.....	99	A11044C.....	95	A1440.....	9	A2202.....	20
51-900-9.50.....	99	A1105.....	1	A1445.....	9	A2204.....	20
51-903-6.25.....	99	A11070.....	95	A1472.....	9	A2206.....	20
51-903-6.50.....	99	A11072.....	95	A1473.....	9	A2208.....	21
51-903-6.75.....	99	A11076.....	95	A1490.....	9	A2210.....	21
51-903-7.00.....	99	A1110.....	2	A1495.....	9	A2211.....	21
51-903-7.25.....	99	A1115.....	2	A1510.....	10	A2238.....	21
51-903-7.50.....	99	A1120.....	2	A1515.....	10	A2240.....	21
51-903-7.75.....	99	A1125.....	2	A1520.....	10	A2242.....	21
51-903-8.00.....	99	A1130.....	2	A1525.....	10	A2244.....	21
51-903-8.25.....	99	A1135.....	2	A1530.....	10	A2246.....	22
51-903-8.50.....	99	A1140.....	3	A1648.....	11	A2247.....	22

Numerical Index

A2248.....	22	A3372.....	29	A3636.....	36	A4073.....	41
A2250.....	22	A3374.....	29	A3638.....	36	A4075.....	41
A2252.....	22	A3376.....	29	A3650.....	36	A4077.....	41
A2254.....	22	A3378.....	29	A3675.....	37	A4079.....	42
A2264.....	22	A3380.....	30	A3677.....	37	A4080.....	42
A2267.....	23	A3382.....	30	A3679.....	37	A4082.....	42
A2268.....	23	A3384.....	30	A3680.....	37	A4084.....	42
A2269.....	23	A3386.....	30	A3681.....	37	A4086.....	42
A2271.....	121	A3388.....	30	A3682.....	37	A4088.....	45
A2273.....	121	A3500.....	31	A3683.....	37	A4112.....	43
A2277.....	23	A3500S	31	A3684.....	37	A4113.....	43
A2278.....	23	A3500T	31	A3686.....	37	A4114.....	43
A2280.....	23	A3502.....	31	A3688.....	37	A4115.....	43
A2286.....	24	A3503.....	31	A3690.....	37	A4160.....	45
A2289.....	24	A3505.....	31	A3692.....	37	A4161.....	44
A2290.....	24	A3506.....	31	A3710.....	101	A4162.....	44
A2294.....	24	A3513.....	31	A3712.....	101	A4163.....	44
A2297.....	24	A3514.....	31	A3714.....	101	A4164.....	44
A2299.....	24	A3515.....	31	A3800.....	63	A4168CX	45
A2300.....	24	A3515T	31	A3801.....	63	A4169.....	45
A2301.....	25	A3516.....	32	A3802.....	63	A4174.....	45
A2302.....	25	A3516T	32	A3804.....	63	A4482.....	27
A2302E	25	A3518.....	32	A3810.....	63	A4482L	27
A2303.....	25	A3520.....	32	A3812.....	63	A4500.....	46
A2304.....	25	A3521.....	32	A3814.....	63	A4502.....	46
A2305.....	25	A3523.....	32	A3820.....	64	A4505.....	46
A2306.....	25	A3526.....	32	A3822.....	64	A4507.....	46
A2307.....	26	A3530.....	33	A3824.....	64	A4510.....	46
A2308.....	26	A3532.....	33	A3830.....	64	A4512.....	46
A2310.....	26	A3534.....	33	A3831.....	64	A4528.....	46
A2310S	26	A3536.....	33	A3832.....	64	A4529.....	46
A2311.....	26	A3537.....	33	A3834.....	64	A4530.....	46
A2312.....	26	A3538.....	33	A3840.....	65	A4531.....	46
A2313.....	26	A3550.....	33	A3842.....	65	A4532.....	46
A2313A.....	26	A3552.....	33	A3844.....	65	A4533.....	46
A2316.....	27	A3554.....	34	A3850.....	38	A4536.....	48
A2316G.....	27	A3556.....	34	A3852.....	38	A4537.....	48
A2316MG.....	27	A3557.....	34	A3860.....	38	A4538.....	48
A2316R	27	A3558.....	34	A3862.....	38	A4539.....	48
A2316RM	27	A3560.....	34	A4000.....	39	A4544.....	47
A2316RS.....	27	A3562.....	34	A4002.....	39	A4545.....	47
A2316T	27	A3564.....	34	A4004.....	39	A4546.....	47
A3320.....	28	A3566.....	34	A4006.....	39	A4547.....	47
A3330.....	28	A3572.....	35	A4011.....	39	A4548.....	47
A3332.....	28	A3574.....	35	A4013.....	39	A4549.....	47
A3334.....	28	A3576.....	35	A4015.....	40	A4550.....	47
A3336.....	28	A3590.....	35	A4017.....	40	A4551.....	47
A3340.....	28	A3592.....	35	A4021.....	40	A4570.....	49
A3342.....	28	A3594.....	35	A4023.....	40	A4572.....	49
A3344.....	28	A3602.....	35	A4025.....	40	A4574.....	49
A3345.....	28	A3624.....	35	A4027.....	40	A4580.....	49
A3346.....	28	A3630.....	36	A4061.....	41	A4582.....	49
A3349.....	28	A3632.....	36	A4062.....	41	A4584.....	49
A3370.....	29	A3634.....	36	A4065.....	41	A4586.....	49

Numerical Index

A4595.....	49	A5080-28.6.....	59	A5804.....	71	A6108.....	80
A4599.....	50	A5081.....	59	A5806.....	71	A6110.....	80
A4600.....	50	A5084.....	59	A5812.....	71	A6162.....	79
A4601.....	50	A5085.....	59	A5828.....	71	A6169.....	79
A4602.....	50	A5086-22	59	A5832.....	72	A6170.....	79
A4604.....	50	A5086-28.6.....	59	A5833.....	72	A6170T	79
A4605.....	50	A5098.....	59	A5834.....	72	A6172.....	79
A4606.....	50	A5100.....	60	A5835.....	72	A6172T	79
A4800CX	54	A5102.....	60	A5854.....	72	A6300.....	81
A4802CX	54	A5110.....	60	A5855.....	72	A6304.....	81
A4804CX	54	A5112.....	60	A5857.....	72	A6312.....	81
A4810CX	54	A5114.....	61	A5884.....	73	A6316.....	81
A4813.....	54	A5116.....	61	A5886.....	73	A6320.....	82
A4814.....	54	A5117.....	61	A5900.....	73	A6322.....	82
A4815.....	54	A5119.....	62	A5908.....	73, 112	A6330.....	82
A4818.....	54	A5120.....	62	A5916.....	73	A6332.....	82
A4819.....	54	A5124.....	62	A5926.....	73	A6342.....	82
A4830.....	55	A5130.....	62	A5927.....	73	A6344.....	83
A4842.....	55	A5132.....	62	A5937.....	74	A6345.....	83
A4900.....	55	A5140.....	62	A5938.....	74	A6346.....	83
A4902.....	55	A5142.....	62	A5940.....	74	A6347.....	83
A4904.....	55	A5400.....	66	A5944.....	74	A6350.....	83
A4906.....	55	A5416.....	66	A5944L	74	A6352.....	83
A4908.....	55	A5450.....	66	A5945.....	74	A6354.....	84
A4949.....	56	A5452.....	66	A5958.....	75	A6356.....	84
A4950.....	56	A5454.....	66	A5960.....	75	A6358.....	84
A4951D	56	A5462.....	66	A5966.....	75	A6370.....	81
A4955.....	56	A5498.....	67	A5967.....	75	A6370P	81
A4959.....	56	A5500.....	67	A5969.....	75	A6376.....	84
A4960.....	56	A5500F	67	A5970.....	75	A6384.....	80
A5002.....	57	A5502.....	67	A5976.....	76	A6388.....	80
A5004.....	57	A5506.....	67	A5978.....	76	A6600.....	85
A5006.....	57	A5510.....	67	A5980.....	76	A6602.....	85
A5040.....	57	A5512.....	67	A5981.....	76	A6604.....	85
A5040P	58	A5542.....	68	A5983CX	76	A6605.....	85
A5042.....	57	A5558.....	68	A5985CX	76	A6640.....	85
A5042P	58	A5559.....	68	A5987CX	77	A6642.....	85
A5044.....	57	A5560.....	68	A5989CX	77	A6644.....	85
A5044P	58	A5562.....	68	A6011.....	77	A6692.....	86
A5050.....	58	A5580.....	69	A6014.....	77	A6693.....	86
A5051.....	102	A5582.....	69	A6017F	77	A6694.....	86
A5052.....	58	A5586.....	69	A6020.....	77	A6816L	99
A5053.....	102	A5588.....	69	A6021.....	77	A6818L	99
A5054.....	102	A5590.....	69	A6022.....	78	A6820L	99
A5055.....	102	A5592.....	69	A6023.....	78	A6822L	99
A5056.....	58	A5656.....	70	A6024.....	78	A6824L	99
A5057.....	102	A5657CX	70	A6025.....	78	A6826.....	99
A5058.....	58	A5660.....	70	A6027.....	78	A6826L	99
A5059.....	102	A5663CX	70	A6028.....	78	A6828.....	99
A5060 - 22	102	A5674.....	70	A6029.....	78	A6828L	99
A5060 - 28.6.....	102	A5674F	70	A6030.....	79	A6832.....	99
A5063.....	102	A5682.....	71	A6034.....	79	A6832L	99
A5064.....	102	A5686.....	71	A6100.....	80	A6834.....	99
A5080-22	59	A5688.....	71	A6102.....	80	A6834L	99

Numerical Index

A6836.....	99	A7033.....	98	A7268.....	104	A7402.....	106
A6836L.....	99	A7034.....	98	A7270.....	104	A7404.....	106
A6838.....	99	A7035.....	98	A7272.....	104	A7406.....	106
A6838L.....	99	A7036.....	98	A7280.....	104	A7410.....	106
A6840.....	99	A7037.....	98	A7282.....	104	A7412.....	106
A6840L.....	99	A7038.....	98	A7286.....	103	A7414.....	105
A6842.....	99	A7039.....	98	A7288.....	103	A7416.....	105
A6842L.....	99	A7040.....	98	A7289.....	103	A7417.....	105
A6844.....	99	A7042.....	98	A7298.....	103	A7418.....	105
A6844L.....	99	A7043.....	98	A7299.....	103	A7419.....	105
A6846.....	99	A7044.....	98	A7300.....	103	A7420.....	105
A6846L.....	99	A7045.....	98	A7302.....	103	A7422.....	105
A6848.....	99	A7046.....	98	A7304.....	103	A7662.....	89
A6848L.....	99	A7047.....	98	A7305.....	103	A7667.....	107
A6850.....	99	A7048.....	98	A7306.....	103	A7669.....	107
A6850L.....	99	A7049.....	98	A7307.....	103	A7669A.....	107
A6852.....	99	A7050.....	98	A7308.....	103	A7671/105.....	107
A6852L.....	99	A7051.....	98	A7309.....	103	A7735.....	107
A6854.....	99	A7052.....	98	A7310.....	103	A7736.....	107
A6854L.....	99	A7053.....	98	A7312.....	103	A7792.....	90
A6856.....	99	A7060.....	88	A7314.....	103	A7799.....	90
A6856L.....	99	A7062.....	88	A7340.....	105	A7804.....	90
A6858.....	99	A7064.....	88	A7342.....	105	A7812.....	90
A6858L.....	99	A7074.....	88	A7344.....	105	A7819.....	90
A6860.....	99	A7077.....	101	A7346.....	105	A7822.....	90
A6860L.....	99	A7078.....	101	A7347.....	106	A7824.....	90
A6862.....	99	A7079.....	101	A7348.....	106	A7826.....	90
A6862L.....	99	A7080.....	101	A7358.....	105	A7827C.....	91
A6868.....	99	A7082.....	101	A7360.....	106	A7827M.....	91
A6868L.....	99	A7102.....	101	A7370.....	105	A7827R.....	91
A6870.....	99	A7104.....	101	A7372.....	105	A7828.....	90
A6870L.....	99	A7106.....	101	A7374.....	105	A7839.....	91
A6872.....	99	A7108.....	101	A7376.....	105	A7840.....	91
A6872L.....	99	A7109.....	101	A7377.....	105	A7841.....	91
A6874.....	99	A7110.....	101	A7378.....	106	A7842.....	91
A6874L.....	99	A7112.....	101	A7379.....	105	A7846.....	92, 110
A6876.....	99	A7122.....	101	A7380.....	106	A7880.....	93
A6876L.....	99	A7124.....	101	A7380L.....	106	A9201.....	94
A6878.....	99	A7126.....	101	A7381.....	106	A9202.....	94
A6878L.....	99	A7130.....	101	A7382.....	106	A9203.....	94
A6880.....	99	A7132.....	101	A7384.....	106	A9210.....	94
A6880L.....	99	A7134.....	101	A7386.....	106	A9220.....	94
A6882.....	99	A7198S.....	101	A7387.....	106	A9222.....	94
A6882L.....	99	A7250.....	104	A7388.....	106	A9230.....	94
A6884.....	99	A7251.....	104	A7388L.....	106	A9232.....	94
A6884L.....	99	A7252.....	104	A7390.....	106	A9239.....	94
A6910.....	87	A7253.....	104	A7391.....	106	A9240.....	94
A6912.....	87	A7254.....	104	A7392.....	105	A9241.....	94
A6914.....	87	A7256.....	104	A7394.....	106	A9242.....	94
A7012S.....	99	A7258.....	104	A7395.....	106	A9246.....	94
A7022S.....	99	A7260.....	104	A7396.....	106	A9248.....	94
A7023S.....	99	A7262.....	104	A7398.....	104	A9304.....	93
A7030.....	98	A7265.....	104	A7400.....	104	A9306.....	93
A7032.....	98	A7266.....	104	A7401.....	104	A9310.....	93

Numerical Index

A9311.....	93	M2658.....	125	M4463.....	138	M5560.....	133
A9312.....	93	M2801.....	123	M4464.....	138	M5561.....	133
A9313.....	93	M2802.....	123	M4467A.....	138	M5562.....	133
A9314.....	93	M2804.....	123	M4471.....	140	M5563.....	133
A9520A.....	108	M2817.....	123	M4474.....	140	M5570.....	134
A9520AC.....	108	M2822.....	123	M4475.....	141	M5571.....	134
AGV-FP7.....	122	M2824.....	123	M4475A.....	141	M5572.....	134
AGV-FP8.....	122	M3000D.....	128	M4477.....	140	M5573.....	134
AGV-FX1.....	122	M3000E.....	128	M4478.....	140	M5580.....	133
AGV-S2.....	122	M3000F.....	128	M4478B.....	140	M5581.....	133
AMD 1020.....	143	M3000G.....	128	M4489.....	142	M5582.....	133
AMD 1021.....	143	M3602.....	127	M4489A.....	141	M5583.....	133
AMD 1061.....	143	M3610.....	127	M4491.....	142	M5584.....	133
AMD 1062.....	143	M3610B.....	127	M4492.....	142	M5590.....	133
AMD 1063.....	143	M3610C.....	127	M4493.....	142	M5591.....	133
C4100.....	121	M3610S.....	127	M4494.....	142	M5592.....	133
C4200.....	121	M3615.....	127	M4495.....	142	M5593.....	133
COX03.....	53	M3626.....	127	M4496.....	142	M5594.....	133
COX04.....	53	M3626A.....	127	M4497.....	142	M6524.....	136
COX05.....	53	M3627.....	127	M4498.....	141	M6575.....	131
COX06.....	53	M3627A.....	127	M4499.....	140	M6575A.....	131
COX07.....	53	M3629.....	127	M4499A.....	140	M6576.....	131
COX2.....	53	M3630.....	127	M4500.....	142	M6576A.....	131
M0275.....	124	M3899.....	130	M4501.....	142	M6594.....	131
M0278.....	124	M3899A.....	130	M4502.....	142	M6595.....	131
M0279.....	124	M3900.....	129	M4503.....	142	M6596.....	131
M0637.....	124	M3902.....	129	M4504.....	141	M6600.....	134
M0638.....	124	M3904.....	130	M4510.....	139	M6602.....	134
M0641.....	124	M3906A.....	130	M4515.....	139	M6605.....	134
M0641A.....	124	M3907.....	130	M4530.....	141	M6611.....	133
M0642.....	124	M3908.....	130	M4531.....	141	M6634.....	133
M0642A.....	124	M3909.....	126	M4619.....	131	M7500.....	123
M0643.....	124	M3910.....	126	M4620.....	131	M7510.....	123
M1000.....	134	M3910A.....	126	M5331.....	136	M8604.....	132
M2270.....	125	M3911.....	130	M5405.....	136	M8604A.....	132
M2270A.....	125	M3911A.....	130	M5501.....	135	M8604B.....	132
M2272.....	126	M3912.....	130	M5502.....	135	M8605.....	132
M2273.....	125	M3914.....	123	M5505.....	135	M8605A.....	132
M2273A.....	125	M3917.....	126	M5510.....	135	M8608.....	132
M2273B.....	125	M3919.....	129	M5511.....	135	M8620C.....	132
M2273C.....	129	M3926A.....	126	M5512.....	135	M8621.....	132
M2273D.....	129	M3927.....	126	M5513.....	135	M8622.....	132
M2273E.....	129	M3992.....	129	M5514.....	135	M8645.....	132
M2273F.....	129	M4410.....	140	M5515.....	135	M9272.....	137
M2273G.....	125	M4447.....	140	M5518.....	135	M9274.....	137
M2273H.....	129	M4448.....	140	M5519.....	135	M9275.....	137
M2274.....	129	M4449.....	141	M5520.....	135	M9276.....	112, 137
M2275.....	125	M4452.....	140	M5521.....	135	M9277.....	137
M2279.....	125	M4453.....	140	M5523.....	135	M9278.....	137
M2280.....	125	M4454.....	140	M5531.....	136	M9280.....	126
M2639.....	125	M4456.....	141	M5535.....	136	M9281.....	126
M2639A.....	126	M4457.....	141	M5536.....	136	M9290.....	100, 137
M2642.....	125	M4458.....	141	M5537.....	135	M9291.....	100, 137
M2657.....	125	M4459.....	141	M5538.....	136	MCD26.....	117

Numerical Index

MCS26.....	117	NST100.....	122
MCU26.....	117	PD-25RFD.....	56
MDP30.....	117	PD.100.00.....	108
MDP35.....	117	SGP06-SG28.....	122
ME-100.....	12	TE.....	122
ME-106.....	12	TI.....	122
ME-120.....	12	V20LCD10.....	90, 108
ME-130.....	12	V28LCD10.....	90, 108
ME-231.....	13	V30PD.....	109
ME-232.....	13	V3MIRD10.....	92, 110
ME-233.....	13	V4MIRD10.....	92, 110
ME-235.....	13	VBCD.....	109
ME-237.....	13	VCAPSD10.....	91, 110
ME-285.....	13	VDGTL1.....	90
ME-420.....	14	VDGTLCF.....	90
ME-425.....	14	VDGTLCM.....	90
ME-451.....	14	VDGTLHM.....	90
ME-701.....	13	VDGTLWF.....	90
MLR50.....	117	VFD.....	109
MLR55.....	117	VG1.....	92
MLR60.....	117	VG3.....	92
MRS32.....	118	VG4.....	92
MRS35.....	118	VG4HAN2.....	92
MRS37.....	118	VHRC.....	91
MRS41.....	118	VHRWF.....	91
MRS52.....	118	VIRIDD10.....	91, 110
MSL**SH.....	116	VMD.....	109
MSL12.....	115	VSLTD10.....	92, 110
MSL14.....	115	VWFD.....	109
MSL16.....	115		
MSL18.....	115		
MSL20.....	115		
MSL22.....	115		
MSL23.....	115		
MSL24.....	115		
MSL25.....	115		
MSL26.....	116		
MSL27.....	116		
MSL28.....	116		
MSL29.....	116		
MSL30.....	116		
MSL32.....	116		
MSL35.....	116		
MST15.....	118		
MST15MK.....	118		
MST22.....	118		
MST30.....	118		
MST45.....	118		
MVR19.....	117		
MVR20.....	117		
MVR20A.....	117		
MVR23.....	117		
MVR23A.....	117		
MVR25.....	117		

Notes

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.